



# OPERATING MANUAL

TA IQ FWD StandUP



This product fulfils the requirements of  
Regulation (EU) 2017/745 on medical devices.

*inspire  
joy of life*

# Contents

<b>Symbols key</b>	6
<b>Introduction</b>	6
<b>Models</b>	6
<b>Indications / contraindications</b>	6
Exit assistance function	7
<b>Delivery and arrival</b>	8
<b>Use only as intended</b>	8
<b>Usage</b>	8
<b>Adjustment</b>	9
<b>Use with products from other manufacturers</b>	9
<b>When the electric wheelchair is used by another person</b>	9
<b>Product life span</b>	10
<b>Base position</b>	10
<b>Overview</b>	11
Model TA IQ FWD StandUp	11
<b>Handling and care of the electric wheelchair</b>	12
Securing the wheelchair	12
Operational inspection	12
Driving characteristics	12
<b>Brakes</b>	12
Service brake	12
Applying the wheelchair's brakes	12
Emergency braking	12
Parking brake	12
Activating the brakes	13
Releasing the brakes	13
Main circuit breakers for all electrical functions	14
Drive/push mode	15
Setting the chair on push mode	15
Setting the chair on drive mode	16

<b>Preparation for use</b>	17
Inspection before driving	19
Charging the battery	19
Position of the control module	21
Description of functions	21
Swing the control module out to the side of the armrest.	21
to drive closer up to a table.	21
Position of the armrest	21
Checking the locking feature	21
<b>Leg supports</b>	22
Footplate	22
Adjusting the footplate height	23
Adjusting the leg support angle	23
<b>Knee pads</b>	24
Removing/inserting the knee pads	24
Adjusting knee pad clearance	25
<b>Armrests</b>	26
Folding up the armrests	26
Folding down the armrests	27
The armrest's locking feature	27
Checking the armrest's locking feature	27
Armrest angle adjustment	28
Adjusting the armrest angle	28
<b>Back support</b>	29
Back cushion	29
<b>Seat</b>	30
Seat cushions	30
Adjusting seating comfort	30
Seat rocker	31
Adjustment of an electronic seat rocker	31
Seat lift	32
Stand up function	33
Reclining function	35
Exit assistance	36

<b>Headrest</b>	38
Adjusting the headrest	38
Use of the headrest in disability vehicles	38
<b>Securing bags</b>	38
<b>Safety belt</b>	39
<b>Safety belt</b>	40
Using the safety belt	40
Opening the safety belt	41
Adjusting belt length	41
<b>Lights</b>	42
<b>Transfer and transport</b>	43
Loading/unloading	43
Ramps and lifts	43
Personal transport in motor vehicles	44
Securing the wheelchair for transport	44
<b>Tyres</b>	45
<b>Maintenance</b>	45
Servicing	45
Service plan	46
Main fuse	48
Lights	49
Running lights	49
Troubleshooting	50
<b>Basic safety precautions</b>	51
Companions	51
Exiting the electric wheelchair	51
Picking up objects	52
Driving on gradients or steep roads	52
Manoeuvring over obstacles	52
Electrical system	53
Carriage in public transport	53
Driving on public roads	54

<b>Cleaning</b>	54
Treated surfaces	55
Disinfection	55
<b>Repairs</b>	56
Assembly	56
Service	56
Replacement parts	56
Precautions for long-term storage	56
<b>Disposal</b>	57
Operating manual for distributors	57
Programming driving characteristics	57
<b>Technical data</b>	58
Air pressure for pneumatic tyres	58
Operating range	58
Driving characteristics when negotiating gradients	58
Applicable standards	59
Specifications in conformance with ISO 7176-15 for model TA IQ FWD StandUp	60
Supplementary technical data for model TA IQ FWD StandUp	62
Explanation of the labels on the electric wheelchair	64
Explanation of the symbols on the nameplate	65
Explanation of the washing instruction symbols	66
<b>DOCUMENTATION OF SERVICING</b>	67
<b>Warranty/guarantee</b>	68
Proof of warranty	69
Service documentation upon sale or transfer	69
<b>Notes</b>	70

## **SYMBOLS KEY**

Safety precautions with a coloured background must always be followed!

-  This symbol designates a precaution or recommendation.
- [ ] Refers to an image number.
- ( ) Refers to a functional component in an image.

## **INTRODUCTION**

Read this manual before you use the wheelchair for the first time and follow the instructions carefully.

If children or teenagers will be using the wheelchair, they should first read this manual along with their parents or guardians before using the chair.

This Operating manual is designed to allow you to enjoy using your electric wheelchair while avoiding accidents.

-  As the wheelchair's accessories may differ from your actual model,

chapters with options that may not be relevant to your particular wheelchair have been included. The order form for your electric wheelchair includes a list of accessories and parts that are available.

PDF files with additional information about our products are available for visually impaired users on our website:

[www.ta-service.dk](http://www.ta-service.dk).

-  Contact your distributor if you have any questions.

Alternatively, visually impaired users can ask a helper to read the text in the documentation aloud.

Stay up-to-date on product safety along with potential recalls of our products on our website:

[www.ta-service.dk](http://www.ta-service.dk).

Our electric wheelchair meets the technical and statutory requirements for medical devices. However, if a serious event nonetheless takes place, please write to us at our e-mail address of 'ta-service@ta-service.dk' and contact the Danish Medicines Agency.

## **MODELS**

These operating instructions apply to the following model:

Model TA IQ FWD StandUP

## **INDICATIONS / CONTRAINDICATIONS**

Seek immediate medical attention if you experience any allergic reactions, redness and/or pressure sores when using this electric wheelchair.

To avoid contact allergies, we recommend only using the wheelchair while wearing clothes.

The numerous features of your electric wheelchair allow it to be used even in case of substantially or severely impaired mobility and ambulation due to structural and/or functional injuries to the lower extremities (including amputation, trauma, musculo-skeletal/neuromusculoskeletal movement disorders, e.g. as a result of:

- paralysis
- loss of limbs (amputation of a leg)
- limb defects/deformations
- joint contractures/joint injuries
- or other conditions

When purchasing a particular product, you should furthermore take into account the disabled individual's physical and mental condition, age, home living conditions and personal circumstances.

Any purchase should always be evaluated and approved on a case-by-case basis by a qualified medical professional (medical device or rehabilitation adviser etc.) and be adapted to the specific circumstances associated with the individual's medical history. Moreover, individuals with any of the issues/conditions described in the contraindications may in certain cases experience mental, emotional or physical difficulties when operating any of the listed models.

The electric wheelchair must not be used in case of:

- Cognitive limitations or intellectual disability that precludes independent use of the wheelchair
- Restrictions in activities of daily living in blind or visually impaired individuals, where visual aids or other assistive devices are unable to compensate for the disability
- Individuals under the influence of medication that has a negative effect on the ability to operate the wheelchair (discuss with your doctor or pharmacist first)
- Conditions that prevent an individual from steering the wheelchair independently
- Severe balance and/or sensory disturbances
- The inability to sit

 Ask your doctor, therapist or authorised distributor about these and any other risks concerning your electric wheelchair.

## Exit assistance function

Seek medical advice immediately for unexpected symptoms that could be associated with using the exit assistance.

The indications/contraindications mentioned in this chapter do not replace the indications/contraindications for use of the electric wheelchair.

We recommend using the exit assistance function for the following indications:

- For people with remaining mobility.

The exit assistance must never be used for:

- The wheelchair user's lack of remaining mobility.

The exit assistance function may only be used with the accompanying person in the case of:

- People who do not have the physical ability to exit the wheelchair without assistance.
- Cognitive limitations or intellectual disability that precludes independent use of the exit assistance.

## **DELIVERY AND ARRIVAL**

All products undergo strict quality control procedures at the factory and are packaged in special cardboard boxes.

- ☞ Nonetheless, we ask that you, immediately upon receipt, inspect your electric wheelchair for any damage that may have occurred during transit – ideally while the hauler is present.
- ☞ Hold onto the wheelchair's packaging in case becomes necessary to return the chair later.

## **USE ONLY AS INTENDED**

Your electric wheelchair is designed to improve your independent mobility both indoors and outdoors.

## **USAGE**

Do not use the wheelchair without the assembled leg supports and armrests!

The wheelchair is exclusively intended to transport a seated individual – it is not designed to withstand or transport other types of tensile or gravitational forces.

The electric wheelchair should only be used on a solid, even surface and can be used in the following manner:

- Indoors (e.g. in flats or outpatient/residential institutions)
- Outdoors (e.g. paved paths in a park)
- Never expose your wheelchair to extreme temperatures or dangerous environmental conditions, such as e.g. strong sunlight, extreme cold or salt water.
- Sand, dust, dirt or other particles may embed themselves in moving parts, causing the components to stop working.

Your electric wheelchair comes with a variety of options that allow you to adjust it to suit your physical requirements.

National regulations may prohibit the wheelchair from being transported on buses, trains, planes or helicopters.

- ☞ Check with the specific carrier to find out about any restrictions.
- ☞ Before your departure, ask your airline about the specific terms and conditions of transport and the baggage requirements that apply to your country of origin and your holiday destination.

The electric wheelchair must be used strictly in accordance with the specifications and limits indicated in the *Technical data* on page 58.

## **ADJUSTMENT**

Ensure that any adjustment, configuration or repair work is only performed by an authorised distributor.

The electric wheelchair can be adjusted to fit individual body sizes. Before you use the wheelchair, your authorised distributor should adjust it and instruct you on its features. This adjustment takes into account both your physical limitations and the environment in which you primarily will be using the chair. Before use, make sure that you are satisfied with how the wheelchair is functioning.

If your distributor revises/updates or substantially modifies your electric wheelchair without using the original manufacturer's parts, this may constitute a resale of the wheelchair. This means that your distributor may need to re-assess the terms and conditions of sale and undertake a new risk analysis.

- ☞ To ensure optimal functioning, we recommend routine inspections of your wheelchair's adjustment as well as a new inspection in the event of any lasting changes in your disorder/disability. An adjustment every 6 months is especially recommended in the case of children and adolescents.
- ☞ We recommend regular medical check-ups so that you can actively use the wheelchair on public transport.

## **USE WITH PRODUCTS FROM OTHER MANUFACTURERS**

In general, any combination of your electric wheelchair with components not supplied by the manufacturer constitutes a material modification of your wheelchair. Contact your authorised distributor to find out whether the manufacturer has approved the combination.

## **WHEN THE ELECTRIC WHEELCHAIR IS USED BY ANOTHER PERSON**

The wheelchair's modular construction allows it to be reconfigured for a new user. Before a new person takes the wheelchair into use, a complete inspection must always be performed.

- ☞ The required hygienic measures, which shall be performed in accordance with a validated hygiene plan, must include disinfection.

The authorised distributor's Service and maintenance manual contains information about the reuse and re-purposing of your wheelchair.

## **PRODUCT LIFE SPAN**

We estimate that the product should have an expected average life span of 5 years insofar as the product is used for the intended purpose and all maintenance and service requirements are adhered to. The product's life span depends on the frequency of use, the environment where it is utilised and the care that it receives. Replacement parts may allow you to prolong its life span. As a rule, replacement parts can still be obtained for up to 5 years after a model has been phased out of production.

- ✎ This stated estimated life span shall not be considered to constitute any expression of the warranty period.

## **BASE POSITION**

When driving downhill, uphill or over obstacles, the settings for seat height, back angle and seat tilt should always be in the base position.

Base position is defined as the following:

- If this is not complied with, there is a risk that the chair flips over.
- Seat tilt in a horizontal position (but no more than 10°).
- Seat height at the lowest setting.

# OVERVIEW

## Model TA IQ FWD StandUp

The overview shows the main components and controls of the electric wheelchair.

No. Name

- (1) Headrest
- (2) Back support
- (3) Armrest
- (4) Seat cushion
- (5) Knee pad
- (6) Calf support
- (7) Footplate
- (8) Anti tip wheel
- (9) Control module
- (10) Front headlight
- (11) Drive wheel
- (12) Drive/push mode lever
- (13) Nameplate
- (14) Swivel castor
- (15) Electronics cover
- (16) Tail light
- (17) Locking armrest handle



# **HANDLING AND CARE OF THE ELECTRIC WHEELCHAIR**

## **Securing the wheelchair**

The electric wheelchair must be secured as follows so that it does not accidentally roll:

1. Move the drive/push mode levers on both sides into the drive setting.
2. Turn off the wheelchair using the control module.

## **Operational inspection**

The electric wheelchair's operations and safety must be checked each time before use.

- 🔊 See the chapter '*Inspection before driving*'.

## **Driving characteristics**

You can control your wheelchair's driving direction and speed via the joystick (up to its configured maximum speed limit).

# **BRAKES**

Brake your wheelchair carefully and in good time. This applies especially when driving toward people or going downhill!

## **Service brake**

The electric motors work as a mechanical brake, stopping the wheelchair smoothly without any jolting until it rolls to a stop.

## **Applying the wheelchair's brakes**

To gently brake the wheelchair, slowly move the joystick backwards to the start position (neutral).

## **Emergency braking**

- 🔊 To stop the wheelchair more suddenly, move the joystick in the opposite direction.

## **Parking brake**

The parking brakes only work when the drive/push mode lever, which controls both motors, is moved to the drive position.

- 🔊 The parking brakes release automatically once the joystick is used.

## Activating the brakes

The wheelchair should not be able to be pushed while the brakes are engaged.

To apply the brakes, move the drive/push mode levers on both sides approximately 45° toward the wheelchair chassis and into the drive (upright) position [1].

- ☞ These levers are designed to be operated by a companion.



## Releasing the brakes

Transfers to/from the electric wheelchair must only be undertaken with the wheelchair turned off and the drive/push mode levers on both sides in the drive setting!

Otherwise, accidentally touching the joystick (power/steering lever) can set the wheelchair in motion. – Risk of an accident!

To release the brakes, move the drive/push mode levers on both sides approximately 45° out into the push position [2].

- ☞ These levers are designed to be operated by a companion.



## Main circuit breakers for all electrical functions

As a standard feature, all TA Service wheelchairs manufactured after 26/05/2021 come equipped with a main circuit breaker for all electrical functions, also known as a Flight kit [3].

Removing the key from the Flight kit allows you to interrupt the wheelchair's power supply [4].

This additional safety feature, which cuts off the power supply, is useful when transporting the wheelchair on a plane or securing it from theft.



## Drive/push mode

The electric wheelchair should only be set on push mode or pushed for the sake of manoeuvring or in emergency situations. This should only be done on flat terrain, with the chair in a standstill position.

- ☞ In push mode, the electromagnetic brakes are deactivated.
  - Thus, the wheelchair can only be braked by switching to drive mode.
- ☞ Grasp the back support at armrest level to manoeuvre the wheelchair.

### Setting the chair on push mode

1. Turn off the wheelchair using the control module. If the wheelchair is turned on while the lever is guided into the push position, an error message appears on the display and an acoustic alarm sounds.
  - ☞ See the '*Control module*' chapter in the Operating manual.
2. Release the brakes [1].
  - ☞ Follow the instructions in *Releasing the brakes* on page 13.

The wheelchair can now be pushed.



## Setting the chair on drive mode

1. Activate the brakes [2].
  - 👉 Follow the instructions in *Activating the brakes* on page 13.
2. Turn on the control module.
  - 👉 See the '*Control module*' chapter in the Operating manual.

Now the electric wheelchair can be driven.



## PREPARATION FOR USE

Perform the following consecutive steps to prepare your electric wheelchair for use.

- ☞ Using the control module, charge the batteries before making your first trip.
1. Complete your preparations for driving.  
Move the levers into drive mode [1] to activate the motors. – Remember to activate the brakes.
- ☞ See *Activating the brakes* on page 13.
2. Check the position of the thermal fuse (2).
    - ☞ The thermal fuse (2) must be depressed.




3. Check the position of the control module.
  - ☞ For ordinary driving, move the control module forward and inward until it settles into the stop position [3].
4. Turn on the control module.
  - ☞ Depress the on/off button (4) on the control module.
  - ☞ See the '*Control module*' chapter in the Operating manual.



## Inspection before driving

Before driving, you should inspect the following:

- ☞ The charge level of the battery
- ☞ The configuration of the maximum speed setting
- See the 'Control module' chapter in the Operating manual.



## Charging the battery

Do not place any objects in the socket aside from the battery charger plug! – Danger of a short circuit!

The batteries must only be charged in a dry environment.

As the battery charger contains hazardous electricity, carefully protect it from warmth, moisture, water, splashing or dripping liquids and impact. – Danger of a short circuit and life-threatening injury!

Make sure that there is adequate space for air to circulate around the charger (do not cover it up) so that any heat generated can dissipate. – Risk of a fire!

Place the battery charger on a solid surface when charging the batteries.

Never place the battery charger on the electric wheelchair seat while charging.

- ☞ When charging, always follow the instructions in the battery charger's operating manual.
1. Park the wheelchair in a safe place.
    - ☞ Follow the instructions in *Securing the wheelchair* on page 12.
  2. Insert the charger plug into the control module socket (1).

3. Turn on the battery charger/plug the battery charger power cord into an appropriate socket.
  - ☞ The charging process has now started.
  - ☞ The battery will not charge if the thermal fuse is turned off (2)!
4. Once the charging process is complete, disconnect the battery charger from the mains and pull the charger plug out of the control module socket.



## Position of the control module

For ordinary driving, move the control module forward and inward [1].

### Description of functions

The push buttons and symbols on the module are all described in detail in the Operating manual in the 'Control module' chapter.

### Swing the control module out to the side of the armrest.

Do not place your hands/arms between the components. – Risk of a crushing injury!

The swing arm [1] allows you to move the control module alongside the armrests.

Enabling you, for example:

to drive closer up to a table.

### Position of the armrest

Risk of an accident if the armrests remain in a vertical or backwards-facing position!

For drive mode, the armrests should be folded down and locked [1].

### Checking the locking feature

Check that the locks are functional by briefly tugging on the armrest.

- 👉 Follow the instructions in *Checking the armrest's locking feature* on page 27.



## LEG SUPPORTS

Before using the leg supports, secure the wheelchair against rolling away unintentionally.

☞ Follow the instructions in *Securing the wheelchair* on page 12.

The angle of the leg support [1] and the height of the footplate can be adjusted via the control module.

☞ See the 'Control module' chapter in the Operating manual for adjustment instructions.

### Footplate

To reduce the overall length, for example to facilitate transport or to make it easier for the user to get in and out of the wheelchair, fold up the footplate [2] [3].

☞ Make sure you don't pinch your fingers!

- Remove both feet from the footplate.

Before driving, the footplate(s) should be folded down again [2].



## Adjusting the footplate height

When adjusting the height of the footplate(s), never touch or handle the adjustment mechanism or the underside of the footplate(s). – Risk of a crushing injury!

Make sure that the leg support and the footplate(s) have adequate ground clearance!

- ⚠ Floor coverings or other surfaces can be damaged when the leg support or footplate(s) are folded down!

Use the control module [1] to raise/lower [1] the leg support to a suitable height.

- ⓘ See the 'Control module' chapter in the Operating manual.

## Adjusting the leg support angle


Never touch or handle the adjustment mechanism when adjusting the angle of the leg support. – Risk of a crushing injury!

Make sure that the leg support has adequate ground clearance!

- ⚠ Floor coverings or other surfaces can be damaged when the leg support or footplate(s) are folded down!

Use the control module to raise [2] or lower [3] the angle of the leg support to a suitable level.

- ⓘ See the 'Control module' chapter in the Operating manual.



# KNEE PADS

Incorrectly adjusted knee pads can cause knee damage in the standing and reclining position.

Check the settings before use and contact your authorised service centre if necessary.

Remove the knee pads before using the exit assistance.


Knee pads that have been removed should not be thrown or dropped, but should be treated properly. – This is the only way to ensure continuous, error-free operation.

The knee pads [1] can be removed and the height and depth adjusted by your service centre.

## Removing/inserting the knee pads

To remove the knee pads, first push in the locking lever (2). Then lift off the knee pads [3] or move the knee pad bracket down the tubing [4].

- ☞ Make sure you don't pinch your fingers!
- ☞ When the locking lever is engaged and released, check the lock.
  - ☞ Pull the knee pads upwards slightly.



## Adjusting knee pad clearance

Do not adjust knee pads positioned by the service centre. – Risk of injury!

Make sure there is one hand width clearance between the leg and the knee pad before you stand up!

☞ Contact your authorised distributor in the event of an unsafe standing position!



To adjust the knee pad clearance, first lift up the locking lever (1).

Then slide the knee pads to the desired distance from the leg [4].

Release the locking lever (1) and allow it to engage in the next possible position to lock the knee pads.

- ☞ Be careful not to pinch your fingers under the raised locking lever!
- ☞ When the locking lever is engaged, check the lock.
  - ☞ To do this, move the knee pads slightly.

# ARMRESTS


Avoid using the armrests [1] to carry or lift the wheelchair.

When the armrest is tilted up, move the support tubing (2) up or down. – Risk of injury!

## Folding up the armrests

To fold up an armrest, depress the button lock (3). This allows the armrest to be guided upwards [4].

- ☞ When the armrest is tilted up, move the armrest support tubing (2) up or down [5]. – Risk of injury!



## Folding down the armrests

For drive mode, the armrests should be swung forward and down [1].

### The armrest's locking feature


Risk of a crushing injury when the armrest engages with the button lock [2]!

1. Swing the armrests forward and down. Move the support tubing downwards and beyond the button lock [2].
2. Press the armrest down onto the button lock [4] until you hear the lock (3) engage.

### Checking the armrest's locking feature

Check that the locks are functional by briefly tugging on the armrest.

- 🔧 Follow the instructions in *The armrest's locking feature* on page 27.
- 🔧 Once you determine that the button lock has engaged, ensure that the locking feature is working properly by tugging on the armrest.



## Armrest angle adjustment

The armrest can be adjusted positively [1] or negatively [2] via a telescopic tube.

### Adjusting the armrest angle

There is a risk of a crushing injury when moving the armrest into the lower position [2]!

1. Pull out the adjusting bolt (3) to adjust the angle of the armrest.
2. Raise/lower the armrest to the desired angle while pulling out the adjusting bolt (3).
3. To lock the armrest into position, push the adjusting bolt (3) into the hole that puts the armrest at the most suitable angle.
4. Move the armrest until you can hear the adjusting bolt (3) engage.

👁️ Check that it is secure by briefly lifting up the armrest!



## BACK SUPPORT

Whenever you adjust the back support, always make sure that the electric wheelchair is on a level surface. Otherwise, there is a risk that the wheelchair may flip over on an incline!

The angle of the back support [1] can be adjusted with the control module.

- ☞ See the '*Control module*' chapter in the Operating manual.

### Back cushion

The back cushion is mounted with a Velcro band to the back support frame and can be removed for cleaning or maintenance [2].

- ☞ More detailed information is available in the separate documentation on back cushions.



# SEAT

## Seat cushions


The seat cushion [1] is mounted to the base-plate with a Velcro band and is removable.

After cleaning or maintenance, put the seat cushion back on and secure it [1]. – Velcro fastening.

## Adjusting seating comfort

The air pressure can be adjusted for improved comfort by opening and closing the valves (2).

- 📖 More detailed information is available in the separate documentation on seat cushions.



## Seat rocker

The seat rocker [3] should only be adjusted while the electric wheelchair is on a flat, horizontal surface. There is a risk that the wheelchair flips over on an incline.

Put your feet on the footplate before adjusting the seat angle and move the leg support to the base position.

Before operating the seat angle, make sure that the surrounding area is free of obstructions. – Risk of injury!

If the seat rocker is adjusted to an angle greater than 14 degrees, the speed of the wheelchair is automatically reduced and lift is limited to approx. 6 cm.

Whenever the back support is adjusted at an angle, there is an increased risk of overturning.

Before driving, ensure that the seat rocker has not been poorly adjusted and that you sit in the chair securely, even on upward or downward sloping terrain.



## Adjustment of an electronic seat rocker

The seat rocker [4] can be controlled either with the control module or a dedicated controller.

- 🔊 Keep your hands and other body parts away from moving parts!
- 🔊 See the '*Control module*' chapter in the Operating manual.

## Seat lift

Make sure that there are no people or obstacles in the area before adjusting the seat height. – Risk of injury!

Do not touch or come into contact with the bottom of the seat before or during seat height adjustment. – Risk of a crushing injury!

Seat height adjustment is only allowed on a flat, level surface.

Do not lower the leg support(s) while the seat is being lifted out of the base position!

The seat lift [1] can be adjusted with the control module.

- ☞ See the 'Control module' chapter in the Operating manual.
- ☞ Once the seat is lifted up from the base position, seat tilt is limited to 14° and at seat lift adjustments greater than 6 cm the speed of the wheelchair is limited!
- ☞ Once the seat returns to the base position, this speed limitation is automatically lifted.



## Stand up function

Make sure that there are no people or obstacles in the area before standing up.  
– Risk of injury!

The stand up function is only allowed when stationary on a flat, level surface.

– The more you stand up, the greater the risk of overturning.

The stand up function must not be used unless a thorough demonstration has been provided by the distributor, therapist or TA Service's product specialist.

An incorrectly adjusted stand up function may result in the user not being able to safely operate the functions of the electric wheelchair.

Both mechanical and electrical adjustments can be made (optimisation of programmable parameters).

Do not lower the leg support(s) while the seat is being lifted out of the base position!



The stand up process [1] can be adjusted with the control module.

- ☞ The user should only stand up with the safety belt [1] and correctly adjusted knee pads [2].
  - ☞ See *Safety belt* on page 40 and *Knee pads* on page 24 for more information.
- ☞ See the '*Control module*' chapter in the Operating manual for adjustment instructions.
- ☞ Make sure you don't pinch your fingers!
- ☞ If the seat has moved up and out of its base position, the speed will be limited.
- ☞ Once the seat returns to the base position, this speed limitation is automatically lifted.



## Reclining function

Before operating the reclining function, make sure that the surrounding area is free of obstructions. – Risk of injury!

The reclining function is only allowed when stationary on a firm, flat, level surface.

Do not lower the foot plate(s) while the leg support(s) are being lifted out of the base position!

The reclining function [1] can be adjusted with the control module.

- ☞ The reclining function should only be used with the safety belt [2] and correctly adjusted knee pads [3].
  - ☞ See *Safety belt* on page 40 and *Knee pads* on page 24 for more information.
- ☞ See the '*Control module*' chapter in the Operating manual for adjustment instructions.
- ☞ Make sure you don't pinch your fingers!



## Exit assistance

Before operating the exit assistance, make sure that the surrounding area is free of obstructions. – Risk of injury!

The exit assistance is only allowed when stationary on a firm, flat, level surface.

The exit assistance function must not be used unless a thorough demonstration has been provided by the distributor, therapist or product specialist.

An incorrectly adjusted exit assistance function may result in the user not being able to safely operate the functions of the electric wheelchair.

Remove the knee pads before using the exit assistance.

The exit assistance [1] can be adjusted with the control module.

☞ See the '*Control module*' chapter in the Operating manual for adjustment instructions.

Before first use, the exit assistance must be adapted and a thorough demonstration provided by the service dealer, therapist or product specialist.

Both mechanical and electrical adjustments can be made (optimisation of programmable parameters).

Your service centre, therapist or product specialist should also familiarise you with the any differing control element functions.



Proceed as follows to exit the wheelchair:

1. Remove the knee pad, if fitted.
  - ☞ Follow the instructions in *Knee pads* on page 24.
2. Perform any adjustments.
3. Turn off the wheelchair.
  - ☞ Follow the instructions in *Securing the wheelchair* on page 12.
- ☞ The following points may require the help of a companion or assistant.
4. Grab something stable to hold on to.
5. Open the safety belt.
  - ☞ Follow the instructions in *Safety belt* on page 40.
6. Exit the wheelchair.

Transfer into the wheelchair in the reverse order.

- ☞ Check that the components required for driving are working correctly.

## HEADREST

The headrest is completely removable and both the height and angle are adjustable.

### Adjusting the headrest


To adjust the height of the headrest or remove it, loosen [2] the tightening handle (1).

### Use of the headrest in disability vehicles

This headrest is approved for transport in motorised disability vehicles!

## SECURING BAGS

A dual hook (3) can be attached to the headrest bracket to hang, for example, a bag.




## ***SAFETY BELT***

Make sure no objects get stuck under the belt! – This prevents unpleasant pressure points.

Aftermarket installation of a seat belt must only be done by an authorised repair centre.

The belt must not be used as a restraint system for securing either the user or the wheelchair during transport in a motor vehicle.



The belt is screwed into place on the outside of the chair with 8 mm bolts (1).

The safety belt is designed to stabilise your sitting posture and prevent you from tumbling forward out of the wheelchair.

When fastening the safety belt, pull both ends forward and fit them into the buckle until you hear an audible click.

To open the safety belt, press the release button and pull the ends of the belt away from one another.

- ⚠ The safety belt can be loosened and must not sit too tight.

## SAFETY BELT

The safety belt must not be used as a restraint system for securing either the user or the wheelchair during transport in a motor vehicle.

Make sure no objects get stuck under the belt! – This prevents unpleasant pressure points.

The optimal position for the safety belt is level with the solar plexus.

The safety belt [1] prevents the user tilting forward, especially when making electrical adjustments such as the stand up function. It also serves as additional stabilisation of the seating position and can be set steplessly according to the user's requirements.

### Using the safety belt

To put on the safety belt, pull the two belts forward and insert the locking halves into each other until you hear a clicking sound (2).

Secure the belt by pulling the lever (3).

- ☞ The safety belt must not be fastened too tightly.
- ☞ Check that the locks are functioning correctly by applying tension!



## Opening the safety belt

To open the safety belt [4], press the release button (2) and pull the ends of the belt away from one another.

## Adjusting belt length

To adjust the belt length, hold the buckle at right angles to the already tightened belt and move it accordingly.




## LIGHTS

If used outdoors or on public roads, the electric wheelchair should be equipped with an LED lighting system (1)+(2).

The lights can be operated via the control module.

- ☞ See the '*Control module*' chapter in the Operating manual.
- ☞ If visibility is poor, especially in the dark, the lights should be turned on so that you can see better and are more easily seen.
- ☞ Make sure that the headlights, turn signals, tail lights and reflectors are not covered by loose clothing or other objects attached to the wheelchair.



## **TRANSFER AND TRANSPORT**

The electric wheelchair must not be lifted by its back support, leg supports, armrest or upholstery!

Keep the seat height and seat tilt in the base position during transport!


Always turn the electric wheelchair off before lifting it!

Any parts that have been dismantled for loading should be stored in a safe place and carefully reassembled before driving.

No special precautions are needed (e.g. fixation at support points) when transporting the detachable parts.

To conserve space, the following measures may be necessary during vehicular transport:

- Fold up the footplate.
- Adjust the back support.



### **Loading/unloading**

The weight of the electric wheelchair is reduced when you remove detachable components.

The wheelchair can be loaded/unloaded by using a ramp or a lift.

### **Ramps and lifts**

Make sure to follow the instructions for the operating manual of the ramp or lift.

Note the manufacturer's specifications for the product.

The ramp's maximum indicated load height must be greater than the height (h) from the loading surface of the particular vehicle to the ground.

The load capacity of the ramp or lift must be greater than the electric wheelchair's allowable total weight.

## Personal transport in motor vehicles

The nameplate on your electric wheelchair indicates whether the chair is approved as a car seat for motor vehicle transport.

- ☞ Learn more about this in *Explanation of the symbols on the nameplate* on page 65.
- ☞ Electric wheelchairs that are not approved as car seats for motor vehicle transport are marked with an additional sticker. – More information is available in *Explanation of the labels on the electric wheelchair* on page 64.

## Securing the wheelchair for transport

The electric wheelchair must now be secured to the anchorage points.

- ☞ Each of the four anchorage points are marked with a symbol [1]+[2].
- ☞ Refer to 'General safety precautions'.
- ☞ When transporting passengers in a motor vehicle, the wheelchair's safety belt system must be used as a restraint device.

According to ISO 7176-19, carriage of passengers is only permissible with the following restraint systems:

- Four-point restraint system in accordance with ISO 10542.

For information on assembly and operation, refer to the manufacturer's operating manual for the particular restraint system.

- ☞ This document accompanies all restraint systems.

Documentation for various restraint systems are also available by clicking the links below the product headings on our website '[www.ta-service.dk](http://www.ta-service.dk)'.



## **TYRES**

The tyres are made of a rubber compound and may leave marks that are difficult or impossible to remove on some surfaces (incl. synthetic, wood or parquet floors, carpets/rugs and floor coverings). We disclaim any liability for any damage to surfaces owing to wear and tear or chemical reactions with the tyres.

## **MAINTENANCE**

Improper or neglectful care and servicing of the vehicle shall limit the manufacturer's liability.

### **Servicing**

The following service plan contains guidelines for servicing the wheelchair.

- 👉 The service plan does not indicate the actual amount of work observed to be needed for the vehicle.

## Service plan

WHEN	WHAT	COMMENT
<b>Before driving</b>	<b>In general</b> Check for proper functioning.	Perform the inspection alone or together with a helper.
	<b>Inspecting the electro-magnetic brake</b> Move the drive/push mode levers on both sides to the 'drive' setting.	Perform the inspection alone or together with a helper.  If the electric wheelchair can be pushed, then the brake is in need of prompt repair by an authorised repair shop. – Risk of an accident!
Especially when driving at night!	<b>Lights</b> Ensure that the lights and reflectors are functioning flawlessly.	Perform the inspection alone or together with a helper.
<b>Every 2 weeks</b> (depending on the overall mileage)	<b>Check the air pressure of the tyres</b> Tyre pressure:  See <i>Technical data</i> on page 58.	Perform the inspection alone or together with a helper.  Use a tyre pressure gauge.
	<b>Adjusting screws, screw sets</b> Ensure that all screws and nuts are tight.	Perform the inspection alone or together with a helper.  Tighten the loose adjustment screws securely.  Visit an authorised repair shop when necessary.
<b>Every 2 months</b> (depending on the overall mileage)	<b>Check the wear on the tyre treads</b> Minimum tread depth = 1 mm	Perform a visual inspection alone or together with a helper.  Contact an authorised repair shop to repair or replace worn or damaged tyres.

WHEN	WHAT	COMMENT
<p><b>Every 6 months</b> (depending on how often the wheelchair is used)</p>	<p><b>Check</b></p> <ul style="list-style-type: none"> <li>- Cleanliness</li> <li>- Overall condition</li> </ul>	<p>See <i>Cleaning</i> on page 54. Perform this alone or together with a helper.</p>
<p>Manufacturer's recommendation: <b>Every 12 months</b> (depending on how often the wheelchair is used)</p>	<p><b>Inspection</b></p> <ul style="list-style-type: none"> <li>- Vehicle</li> <li>- Battery charger</li> </ul>	<p>Performed by the distributor.</p>


## Main fuse

The thermal fuse's fail-safe button must be pushed in.

The main fuse consists of a thermal fuse (1) that pops out in the event of a current overload.

If the fuse pops out, it must be pushed in again.

- 👉 If the fuse continues to trip, take the vehicle to an authorised service centre to assess the cause of the issue.
- 👉 See *Technical data* on page 58.



## Lights


The bulbs (1) + (2) are long-life LED bulbs.

- ✎ Any defective LED bulb must be repaired by an authorised repair shop immediately.

## Running lights

The lights must be adjusted so that the lower edge of the beam lies approx. 3 metres ahead of the wheelchair with the seat tilt down in the base position.

- ✎ Visit an authorised repair shop for adjustment whenever necessary.



## Troubleshooting

Malfunction	Cause	Solution
The control module's battery indicator does not light up when switched on.	The main thermal fuse has popped out.	Press the thermal fuse back in. If it pops out again, contact your service workshop for repair.
	The power plug is not connecting with the power supply.	Inspect the connectors.
When the wheelchair is turned on, the battery indicator blinks and sounds an acoustic alarm.	One or both of the drive motors are set on push mode.	Move the drive/push mode levers on both sides to drive position.
	The connector on of the drive motors is not properly connected.	Inspect the connectors.
	Electronics issue.	Must be repaired at an authorised repair shop. Push mode. Drive/push mode lever in drive position.
	Other malfunctions.	Refer to 'Troubleshooting' in the control module operating manual.
Lights are not working.	Defective LED bulb.	Let an authorised repair shop perform the repair/replacement.
	Defective control or lighting module.	Let an authorised repair shop perform the repair/replacement.

## **BASIC SAFETY PRECAUTIONS**

These safety precautions are an excerpt of the *Safety precautions and general instructions for use*, which is available on our website: [www.ta-service.dk](http://www.ta-service.dk).

Do not stick your fingers inside any open wheelchair frame tubing (e.g. when leg supports are removed or the armrest is tilted up). – Risk of injury!

When using the electric wheelchair, especially when stopped or on an ascending/descending slope, assume a safe sitting posture. – Risk of an accident!

To assume a safe sitting posture, your back should be flush with the back cushion and your pelvis should be situated all the way back on the seat.

Transfers from the wheelchair should only be done on sloping terrain in emergencies, with the help of a companion and/or assistant! – Risk of an accident!

The tilt should only be adjusted while the electric wheelchair is on a flat, horizontal surface. Otherwise, there is a risk that the wheelchair may flip over on an incline!

If the angle of the back support is adjusted, the risk of overturning increases.

Before driving, make sure that the seat tilt has not been set improperly and that it facilitates a safe sitting posture.

Before and during any trip, ensure that the leg support(s) have adequate ground clearance. – Risk of an accident!

You should not smoke while using the electric wheelchair.

In direct sunlight, the seat cover/cushion, armrest cushions, leg supports and handles/levers can reach temperatures above 41 °C. – This may injure unprotected skin! Park in the shade to prevent the surfaces on your electric wheelchair from getting hot.

The bag holder (accessory), which can be attached to the headrest bracket, consists of special hooks that allows you to hang e.g. a bag. – The maximum tolerable overall load on the bag hooks is 5 kg.

Only enter/exit the wheelchair when it is turned off and the drive/push mode lever is in drive position!

– Otherwise, accidentally touching the joystick can set the wheelchair in motion! – Risk of an accident!

### **Companions**

Your companion should be made aware of any hazardous situations before helping you. Any components that your companion handles on the wheelchair must be securely fastened.

### **Exiting the electric wheelchair**

Drive the wheelchair as close as possible to where you plan to exit the chair.

☞ In connection with this, the precautions in *Securing the wheelchair* on page 12, *Footplate* on page 22 and *Folding up the armrests* on page 26 must be followed.

☞ We recommend that you have someone assist you when exiting the chair.

## Picking up objects

Avoid bending your upper body far forward, to the side or backwards, especially when picking up or setting down heavy objects.

- Risk of the electric wheelchair overturning, especially with narrow seat widths and high seat heights (seat cushion)!

## Driving on gradients or steep roads

As the gradient of the road increases past a certain point, the stability of the wheelchair, braking and steering become increasingly affected due to reduced traction, and overturning becomes a significant hazard.

 See *Technical data* on page 58.

Never lean into the lower side of slopes and road surfaces or paths with shoulder gradients.

Always drive up and down hills at a low speed.

Drive at an appropriate speed when going up- or downhill.

Avoid switching to push mode when ascending or descending a slope. The automatic brakes do not work when the wheelchair is in push mode.

Avoid pushing the wheelchair up or down a slope.

There is a risk of overturning when driving through curves or turning on an upward or downward slope.

Avoid driving on upward or downward gradients that are poorly paved. Even a small amount of, among other things, ice, water, moss or algae on the road may cause the wheelchair to lose its traction and slide uncontrollably. In such cases, put the joystick in neutral immediately.

The braking force that can be applied to the road is substantially lower when driving down a hill compared to driving on a level road and deteriorates further if the road conditions are poorer (e.g. wet surfaces, snow, gravel, sludge). Apply the brakes carefully and in a controlled manner to avoid hazardous skidding and loss of control of the wheelchair.

When approaching the bottom of the hill, to prevent a hazardous sudden stop, make sure the footplate(s) are high enough to clear the ground.

Road surfaces with a side gradient (e.g. pavements with a gradient on each side) cause the wheelchair to drift toward the lower side of the surface. You or your companion can compensate for this by steering in the opposite direction.

## Manoeuvring over obstacles

The ability of the wheelchair to manoeuvre over obstacles is correlated, among other things, with the slope of the road surface and the position of the leg support.

Manoeuvring over an obstacle is hazardous!  
– Danger of overturning the wheelchair!

Manoeuvring over an obstacle is a special risk situation, where you must take into account a number of different safety precautions involving e.g. slopes and driving surfaces with side gradients.

Drive around obstacles, such as cracks or grooves in the roadway, railway tracks, sewer covers or other road hazards, as much as possible.

Always drive toward smaller obstacles, such as ledges/edges, slowly and at a straight angle (90°). Drive forward and approach to within approx. 0.5 m of the obstacle so that the front and rear tyres pass over the obstruction at the same time. If you fail to do this, your wheelchair may fall over and eject you!

Maintain a safe enough distance between yourself and drop-offs, stairs, slopes, embankments etc. that you have sufficient space to react, brake and turn.

If possible, get one or more helpers to lift you out of the wheelchair and carry you to your destination.

If you drive down a step (e.g. a pavement curb), it's easy to fall out of the wheelchair if your footplate(s) and leg support hit the road surface. Be extremely attentive when crossing railway tracks. – Accidental deviation from your course!

It is impossible to navigate stairs safely with ordinary electric wheelchairs.

## Electrical system

Improper and/or inappropriate modifications to the driving characteristics may affect the wheelchair and its safety. – Risk of an accident!

The electrical control system of the wheelchair must never be tampered with.

If the electric wheelchair is acting strangely or driving uncontrollably, put the joystick in neutral immediately and/or turn off the wheelchair.

## Carriage in public transport

Your electric wheelchair is not designed for passenger transit in public transit vehicles. There may be limitations. We recommend using a seat that is affixed to the public transit vehicle.

If seated transport in the wheelchair is nonetheless unavoidable, the following precautions must be observed:

- Use the space reserved by the public transit company for disabled passengers.
  - Before parking the wheelchair, ensure that you are compliant with the regulations of the transit company.
  - Park your wheelchair in the reserved space so that it is facing the direction of travel.
  - The wheelchair should be parked in the reserved space with the backrest flush against a supportive structure.
  - One side of the wheelchair should be likewise be supported by an immovable structure so that the chair does not slide during an accident or sudden braking manoeuvre.
  - Ensure that the motors are set on drive and that the parking brake is engaged.
-  Follow the instructions in *Parking brake* on page 12.

## Driving on public roads

Comply with all national regulations for driving on public roads and ask your authorised distributor about any required accessories.

The electric wheelchair can be delivered with accessory lights. The system consists of:

- lights
- retroreflectors

In poor visibility or at night, we recommend installing the lights and driving with them on, both to improve your visibility and make you more easily seen.

- ☞ When driving on public roads, the driver is responsible for the functional and operational safety of the wheelchair.
- ☞ Always comply with any applicable traffic regulations when driving on public roads.
- ☞ If you drive at night, we recommend wearing bright, conspicuous clothing so that you are more easily seen.
- ☞ Avoid using carriageways and bike paths when driving at night.
- ☞ Make sure that the lights are not covered by loose clothing or other items attached to the wheelchair.

## CLEANING

Non-ionic detergents, solvents and especially alcohol may react with synthetic coverings.

Never rinse or pressure rinse your electric wheelchair. – Danger of a short circuit!

Cushions and coverings are usually supplied with cleaning instructions (care tag).

☞ See *Explanation of the washing instruction symbols* on page 66.

Otherwise, the following precautions apply:

- ☞ Clean the cushions with warm water and washing-up liquid.
- ☞ Stains can be removed with a sponge or soft brush.
- ☞ Use ordinary detergent to remove stubborn stains.
- ☞ If the wheelchair is equipped with a Varilite Back System, both the seat cover and back cushion can be machine washed. Follow the instructions.

Rinse with clean water and leave to dry.

Clean the chassis and the wheels with a mild cleaning agent, then dry thoroughly.

- ☞ Check the frame for corrosion and other damage.
- ☞ Synthetic materials must only be cleaned with warm water and a neutral detergent or brown soap.
- ☞ If you use an ordinary synthetic detergent, follow the manufacturer's instructions.

Always keep the lights clean and check that they are functioning properly before driving.

- ☞ Always protect the electrical components from water and moisture!
  - Sprayed water can damage the electronics and control module.

Water-soluble, silicon-free detergents and cleaning products should be used.

- ☞ Follow the instructions from the product's manufacturer.

Do not use aggressive cleaning agents, such as solvents, or stiff brushes etc.

More detailed information about cleaning and care is available on our website:

*'[www.ta-service.dk](http://www.ta-service.dk)'.*

## Treated surfaces

The unique surface coating on the wheelchair provides the best possible protection against corrosion.

Occasional lubrication of the moving parts with a small amount of lubricating oil, can improve their performance.

## Disinfection

If the product is used by multiple individuals (e.g. at a care institution), ordinary disinfectant should be used.

- ☞ Clean the cushions and handles/levers prior to disinfecting them.
- ☞ Spraying and wiping with an approved, recognised disinfectant is permitted.

Information about tested, recognised disinfectants and disinfection procedures can be obtained from your national authority for public health protection.

- ☞ Disinfectants may damage certain surfaces and impair the performance of components over time.
- ☞ Follow the instructions from the product's manufacturer.

# REPAIRS

In principle, repairs should only be performed by authorised distributors.

## Assembly

You can have complete confidence in your distributor when it comes to repairs. He or she has been instructed on how to perform this work.

## Service

If you have any questions or need assistance, please contact your distributor who can advise you and carry out service and repairs.

## Replacement parts

Replacement parts can only be obtained from your distributor. In case of repair, only original replacement parts may be used!

☞ Components from other manufacturers may cause malfunctions.

Your distributor has a list of replacement parts with associated part numbers and schematics.

To ensure you receive the right replacement parts, always provide the correct serial number (SN) of your electric wheelchair! This number can be found on the nameplate.

For all repairs performed on the wheelchair by the service centre, supplemental information, such as assembly/operating information from the Operating manual, must be provided, and the date of any modifications must be noted and passed on when replacement parts are ordered.


This prevents inaccurate orders when replacement parts are subsequently ordered.

## Precautions for long-term storage

The following precautions are necessary if the wheelchair will not be used for a longer period of time:

- ☞ Charge the batteries for 16 hours at least once a month.
- ☞ The prescribed storage temperatures must be ensured.
- ☞ Follow the instructions in *Technical data* on page 58.

## DISPOSAL



The wheelchair must be disposed of in accordance with applicable national regulations.

Contact your municipality for information about your local environmental recycling centre.

### Operating manual for distributors

On our website '[www.ta-service.dk](http://www.ta-service.dk)' you can find a Service and maintenance manual for this electric wheelchair that includes the following information:

1. Configuration that can be performed with tools
2. Step-by-step descriptions of procedures for important repairs
3. Notes about changes to specific models
4. A checklist for the annual inspection

This includes mandatory functional checks and guidelines associated with the performance of the inspection work.

☞ However, the list does not cover the actual care required by the vehicle.

Once the annual inspection has been properly performed, document the completed inspection in the Operating manual.

A template for additional inspection documentation can be copied from the Service and maintenance manual as needed. If this is utilised, it should be attached to the Operating manual.

### Programming driving characteristics

The electric wheelchair's driving characteristics can be configured via a programming unit.

☞ Please refer to the '*Service and maintenance manual*'.

The wheelchair's driving characteristics should routinely be adapted to the particular user's needs and skill level.

☞ The programming must be customised for the user. The individual's reflexes, general health condition and physical and mental capabilities must all be taken into account. Consultation with the user's doctor or therapist can be extremely helpful in this context.

☞ Any modification from the manufacturer's default programming settings may increase the risk of accidents.

☞ Potential risk of overturning!

## TECHNICAL DATA

All the information in the 'Technical data' chapter applies to the standard model.

Measurement tolerances  $\pm 15$  mm,  $\pm 2^\circ$

Calculating the max. allowable user weight:

The overall load capacity can be calculated by adding the weight of the empty electric wheelchair to the max. allowable user weight.

Added weight due to aftermarket parts or baggage lowers the max. allowable user weight.

Example:

A user wants to take along baggage weighing 5 kg. This reduces the max. allowable user weight by 5 kg.

### Air pressure for pneumatic tyres

The max. air pressure is listed on both sides of the tyre.

Air pressure - swivel castor

Standard:

2.0 bars = 29 psi

Air pressure - drive wheel

Standard:

2.5–3.0 bars = 33–44 psi

### Operating range

Our indicated nominal data are realistic, provided that ISO 7176-4 is fully complied with.

The operating range largely depends on the following factors:

- Condition of the battery
- User weight
- Driving speed
- Driving behaviour
- Road condition
- Driving conditions
- Ambient temperature

The operating range is severely limited by:

- Frequent driving up ramps
- Batteries with deficient charging capacity
- Low ambient temperatures (e.g. in winter)
- Frequent starting and stopping (e.g. in shopping centres)
- Ageing sulphated batteries
- Sudden forced steering manoeuvres
- Reduced driving speed (especially at pedestrian speeds)

In practice, the operating range that usually can be achieved under normal conditions is typically reduced by 40–80% relative to the nominal value.

### Driving characteristics when negotiating gradients

For safety reasons, the wheelchair should only be operated without a driver on upward or downward slopes exceeding the allowable limits (e.g. on ramps)!

## Applicable standards

The electric wheelchair fulfils the following standards:

- EN 12184: 2014
- ISO 7176-8: 2014
- ISO 7176-19: 2008
  - ☞ Crash test evaluations, where the wheelchair was attached via the vehicle's braking system, were performed in accordance with the testing methods in Annex D.
- Four-point restraint system in accordance with ISO 10542.
  - ☞ The applicable product documentation accompanies the products.

The models are categorised as Type class B pursuant to the EN 12184 standard.

Our modules and components fulfil the requirements of EN 1021-2 for ignitability resistance.

## Specifications in conformance with ISO 7176-15 for model TA IQ FWD StandUp

	min.	max.
Overall length (measured with a 0° seat tilt)	985 mm	1160 mm
Overall width	630 mm	720 mm
Total weight, max. allowable		320 kg
User weight (incl. payload) for reclining function		140 kg 100 kg
User weight When the product is used for seating in a motor vehicle (crash tested in accordance with ISO 7176-19)		136 kg
Weight of the heaviest part		26 kg
Actual seat depth	250 mm	590 mm
Actual seat width	370 mm	550 mm
Seat surface height at front edge (without seat cushion) with 0° seat tilt	440 mm	740 mm
Electronic seat angle	0°	45°
Electronic seat lift		300 mm
Angle of electronic back support (measured based on a vertical position from the seat plate)	90°	180°
Back support height	540 mm	665 mm
Thigh control for seat (Thigh length, measured without the seat cushion)	370 mm	580 mm
Static stability driving downhill	15.5°	19.6°
Static stability driving uphill	19.6°	19.6°
Static stability side-to-side	13.5°	19.6°
Dynamic stability driving uphill		10°
Angle - leg supports/seat surface	90°	180°
Height of armrest from seating surface (without seat cushion)	185 mm	285 mm
Back support to front edge of armrest	370 mm	475 mm
Obstacle clearance		100 mm
Min. turning radius (measured with 0° seat tilt)	650 mm	
Weight of test dummy (ISO 7176-8)		140 kg
Top forward driving speed (depends on equipment)	6 km/h	12 km/h

	min.	max.
Minimum braking distance at top speed	2620 mm	
Operating range (at 6 km/h) <small>(depends on battery capacity)</small>		40 km
Operating range (at 10 km/h and 12 km/h) <small>(depends on battery capacity)</small>		35 km
Axle, horizontal position	- mm	- mm

## Supplementary technical data for model TA IQ FWD StandUp

	min.	max.
Noise level		70 dB(A)
Enclosure class		IPX4
Min. turning radius	1170 mm	
Drive controller output		24 V/120 A
Engine output		2 x 350 W
Main fuse		80 A
Lights (accessory)		24V LED type
Payload		5 kg
Front axle pressure (max. allowable)		210 kg
Rear axle pressure (max. allowable)		150 kg
Ground clearance		70 mm
Empty weight (with batteries)		180 kg
Empty weight (without batteries)		129 kg
Overall height	1000 mm	1100 mm

### Transport dimensions

Length (footplates folded up)	820 mm	
Width	630 mm	720 mm
Height <small>(Back folded onto the seat and seat cushion removed from the seat plate and placed on top of the back support)</small>	700 mm	

### Thermodynamic data

Environmental temperature range	-20 °C to +50 °C
Storage temperature with batteries	-20 °C to +50 °C

### Swivel castor

200 x 50 mm diameter (8")	Tyres, max. 2.0 bars (29 psi/200 kPa) puncture-proof
---------------------------	---

	min.	max.
<u>Drive wheel</u>		
364 x 75 mm diameter (14" x 3.5")	Tyres, max. 2.5–3.0 bars (33 -44 psi/250 kPa) puncture-proof	


<u>Batteries</u>		
2 x 12 V 75 Ah (5 h) / 80Ah (20 h)	Enclosed, maintenance-free	
Max. battery dimensions (L x W x H)	260 x 168 x 215 mm	
Charging current		8 A

## Explanation of the labels on the electric wheelchair



### Attention!

Read the Operating manual and attached documentation.



Do not lift the electric wheelchair by the armrests or leg supports.

Do not lift the wheelchair by any detachable parts.



Drive mode



Push mode



Only push the wheelchair on a level surface.



Information about the charging plug



This electric wheelchair is **not** approved for use as a motor vehicle seat.



Warning about the risk of a crushing injury – Do not place your hands/arms between the components.



Operation in public transit vehicles – Turn off or switch to driving program 1.



Max. allowable user weight if approved for use as a motor vehicle seat.

## Explanation of the symbols on the nameplate



Manufacturer



Order number



Serial number



Date of production



Allowable user weight



Allowable total weight



Allowable axle pressure



Allowable upward gradient



Allowable downward gradient

max. ... km/h

Max. allowable speed



The product is approved for use as a car seat.



Max. allowable user weight if approved for use as a motor vehicle seat.



The product is **not** approved for use as a car seat.



Medical device

## Explanation of the washing instruction symbols

(The symbols are in compliance with European standards.)



Wash on gentle cycle at the maximum indicated temperature (°C).



Wash on normal cycle at the maximum indicated temperature (°C).



Hand wash only.



Do not bleach.



Not suited for the dryer.



Do not iron.



Do not dry clean.

# DOCUMENTATION OF SERVICING

## Vehicular data:

Model:

Packing slip no.:

Serial no. (SN):

## Recommended safety inspection: 1st year (after no later than 12 months)

Stamp of speciality dealer:

Signature: \_\_\_\_\_

Location, date: \_\_\_\_\_

Next safety inspection in 12 months

Date: \_\_\_\_\_

## Recommended safety inspection: 2nd year (after no later than 12 months)

Stamp of speciality dealer:

Signature: \_\_\_\_\_

Location, date: \_\_\_\_\_

Next safety inspection in 12 months

Date: \_\_\_\_\_

## Recommended safety inspection: 3rd year (after no later than 12 months)

Stamp of speciality dealer:

Signature: \_\_\_\_\_

Location, date: \_\_\_\_\_

Next safety inspection in 12 months

Date: \_\_\_\_\_

## Recommended safety inspection: 4th year (after no later than 12 months)

Stamp of speciality dealer:

Signature: \_\_\_\_\_

Location, date: \_\_\_\_\_

Next safety inspection in 12 months

Date: \_\_\_\_\_

## Recommended safety inspection: 5th year (after no later than 12 months)

Stamp of speciality dealer:

Signature: \_\_\_\_\_

Location, date: \_\_\_\_\_

Next safety inspection in 12 months

Date: \_\_\_\_\_

## **WARRANTY/GUARANTEE**

Any lack of compliance with the Operating manual or unprofessionally performed maintenance work, especially with respect to technical modifications and add-ons (installed parts) without our express consent shall not only void this warranty/guarantee, but also nullify our general product liability.

The product comes with a 2-year warranty.

However, the warranties for the batteries and charger are 1 year each.

The warranty period begins on the date of purchase.

For any warranty/guarantee claims, please contact your distributor and present the following PROOF OF WARRANTY certificate along with any necessary information about the model, packing slip number, date of delivery and serial number (SN).

The serial number (SN) is visible on the nameplate.

Acceptance of warranty/guarantee claims is contingent on the proper use of the product, the use of original replacement parts from an authorised distributor and the routine performance of maintenance and inspection.

The warranty/guarantee does not cover surface damage, wheel inner tubes or tyres, damage due to loose nuts or screws, or damaged bores and screw holes resulting from repeated assembly of the chair.

Likewise, the warranty/guarantee does not cover damage to motors or electronics due to unprofessional cleaning with steam-cleaning equipment or deliberate or accidental exposure of the components to water.

No liability shall be assumed for defects or malfunctions due to sources of radiation such as high-transmission mobile phones, stereo equipment and other powerful sources of noise or interference in excess of standard specifications.

This Operating manual must be passed on to any new owner or user of the product.

Please refer to our website '[www.ta-service.dk](http://www.ta-service.dk)' to evaluate our products.

We reserve the right to make technical modifications for the purpose of product improvement.

## Proof of warranty

Please complete! Can copy this as needed and send it to your distributor.

# Warranty/guarantee

Model name:

Packing slip no.:

SN (refer to the nameplate):

Date of delivery:

Distributor's stamp:

## Service documentation upon sale or transfer

### Vehicular data:

Serial no. (SN):

Model:

Packing slip no.:

Stamp of speciality dealer:

Signature: \_\_\_\_\_

Location, date: \_\_\_\_\_

Next safety inspection in 12 months

Date: \_\_\_\_\_

# **NOTES**

---





### **Manufacturer:**

📍 TA Service A/S  
Centervej Syd 2  
DK-4733 Tappernøje  
Denmark

✉ ta-service@ta-service.dk  
T: +45 56 72 57 77  
www.ta-service.dk

---

### **Distributor**

Please find local distributors at  
[www.ta-service.dk](http://www.ta-service.dk)