

## TA Indoor Wave


<b>COMPANY:</b>		<b>Shipping address:</b>	
Customer number:	Date:	Address:	
		Postcode:	Town/City:
Cost Estimate:	Order:	Country:	Phone:
Contact person:		E-Mail:	
Reference:			

### Specifications

Lenght with leg support: (mm)	885
Overall width: (mm)	560
Seat width min.-max.: (mm)	250-550
Seat depth min.-max.: (mm)	250-600
Seat height w/o cushion: (mm)	400
Backrest height min.-max.: (mm)	540-665
Seat tilt: (°)	45
Backrest angle min.-max.: (°)	50-96
Lift: (mm)	320
Footplate to seat distance min.-max.: (mm)	370 - 580
Leg to seat surface angle min.-max.: (°)	90-180
Arm support to seat distance min.-max.: (mm)	185 - 285

Total mass: (kg)	113,5
Maximum user weight: (kg)	125
Dynamic stability uphill: (°)	10
Obstacle climbing: (mm)	65
Ground clearance: (mm)	700
Maximum speed forward: (km/h)	€ 6,00
Range w. batteries 38Ah: (km)	22
Engines:	2pcs. 24 V / 450 Watt
Minimum turning radius: (mm)	450

### User details



All seat dimensions are always measured without seat cushion!

01 Seat depth:
02 Seat width:
03 Chest width:
04 Chest depth:
05 Elbow height:
06 Arm pit height:
07 Back height:
08 Head height:
09 Forearm length:
10 Lower leg length:
Height:
Weight (kg):

Specify dimensions in mm

Refer to the sketches!

The grey colored rows are indicative and show a standard configuration

## 1 Chassis / Motor

CODE	DESCRIPTION		Speed	Motor Power	Range	Max. User weight
✓ 95346	TA Indoor Wave Multifunctional frontwheel electric wheelchair	Lift 320 mm & tilt 45° included	6 km/h	2 x 450 W	15 km	125 kg



95346

## Options for chassis / Motor

CODE	DESCRIPTION
✓ 95025	Cover black
<input type="checkbox"/> 95027	Cover titanium
<input type="checkbox"/> 95350	TA 4 point car attachment
<input type="checkbox"/> 95351	Dahl car docking with lock plate
<input type="checkbox"/> 95352	TA 4 point & Dahl car attachment



95025



95027



95350



95351

## 2 Battery and chargers

CODE	DESCRIPTION	
✓ 95542	Battery, Celpower, 38A IW (2 pcs)	€ 164,00
✓ 95387	Charger 24V / 6A	€ 135,00

## 3 Wheels



### Drive wheels

260 x 85

Rims alu gray

Pneumatic tyres		Puncture-proof tyres	
Tire gray	Tire black	Tire gray	Tire black
<input type="checkbox"/> 95266	<input type="checkbox"/> 95268	<input type="checkbox"/> 95269	<input type="checkbox"/> 95267

### Castor wheels

200 x 50

Rims alu gray

Pneumatic tyres		Puncture-proof tyres	
Tire gray	Tire black	Tire gray	Tire black
<input type="checkbox"/> 97292	<input type="checkbox"/> 97260	<input type="checkbox"/> 97227	<input type="checkbox"/> 97261
<input type="checkbox"/> 97245	<input type="checkbox"/> 97259	<input type="checkbox"/> 97262	<input type="checkbox"/> 97195

Rims alu black

## 4 Seat

Lift / tilt

CODE	DESCRIPTION
✓ 95011	Acutator lift 320 mm
✓ 95546	Seat tilt 45°



95011



95546

### Seat width

CODE	DESCRIPTION
<input type="checkbox"/>	Seat width 250mm adjustable to 250 - 300 or 350 mm (80291-2)
<input type="checkbox"/>	Seat width 380mm including depth extender (80233+80291-4)
<input type="checkbox"/>	Seat width 400mm including depth extender (80233+80291-4)
<input type="checkbox"/>	Seat width 450mm including depth extender (80216+80291-4)
<input type="checkbox"/>	Seat width 500mm including depth extender (80215+80291-4) reinforced with middle plate (80502)
<input type="checkbox"/>	Seat width 550mm including depth extender (80215+80291-4) reinforced with middle plate (80502)

### Seat depth

Seat depth is stepless adjustable from 250 - 600mm:

Specify required seat depth:

-----

### Seat system



80519

CODE	DESCRIPTION	Seat width				
		380mm	400mm	450mm	500mm	550mm
<input type="checkbox"/>	80519 TA Strap seat	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

## Seat cushion



Varilite Evolution™  
PSV

For wheelchair users with a high degree of tissue breakdown and symmetric positioning needs.

VARILITE's PSV takes guesswork out of cushion adjustment and makes it automatic. Evolution PSV uses VARILITE's revolutionary new Pressure Setting Valve (patent pending). The PSV works by sensing the amount of internal cushion pressure. As air is released, internal pressure decreases and the user is immersed in the cushion for effective pressure distribution.

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)
<input type="checkbox"/>	93000	Varilite Evolution™ PSV	300x350
<input type="checkbox"/>	93001	Varilite Evolution™ PSV	350x350
<input type="checkbox"/>	93002	Varilite Evolution™ PSV	350x400
<input type="checkbox"/>	93420	Varilite Evolution™ PSV	350x450
<input type="checkbox"/>	93421	Varilite Evolution™ PSV	350x500
<input type="checkbox"/>	93422	Varilite Evolution™ PSV	380x380
<input type="checkbox"/>	93423	Varilite Evolution™ PSV	380x430
<input type="checkbox"/>	93424	Varilite Evolution™ PSV	400x350
<input type="checkbox"/>	93003	Varilite Evolution™ PSV	400x400
<input type="checkbox"/>	93004	Varilite Evolution™ PSV	400x450
<input type="checkbox"/>	93005	Varilite Evolution™ PSV	400x500
<input type="checkbox"/>	93006	Varilite Evolution™ PSV	430x430
<input type="checkbox"/>	93007	Varilite Evolution™ PSV	450x400
<input type="checkbox"/>	93008	Varilite Evolution™ PSV	450x450
<input type="checkbox"/>	93009	Varilite Evolution™ PSV	450x500
<input type="checkbox"/>	93435	Varilite Evolution™ PSV	500x400
<input type="checkbox"/>	93268	Varilite Evolution™ PSV	500x450
<input type="checkbox"/>	93010	Varilite Evolution™ PSV	500x500
<input type="checkbox"/>	93277	Varilite Evolution™ PSV	550x450
<input type="checkbox"/>	93278	Varilite Evolution™ PSV	550x500
<input type="checkbox"/>	93279	Varilite Evolution™ PSV	600x450
<input type="checkbox"/>	93280	Varilite Evolution™ PSV	600x500



Varilite Meridian™  
PSV

Dual-chambered, adjustable skin protection cushion. Ideal for prevention of sacral sitting

VARILITE's PSV takes guesswork out of cushion adjustment and makes it automatic. Evolution PSV uses VARILITE's revolutionary new Pressure Setting Valve (patent pending). The PSV works by sensing the amount of internal cushion pressure. As air is released, internal pressure decreases and the user is immersed in the cushion for effective pressure distribution.

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)
<input type="checkbox"/>	93362	Varilite Meridian™ PSV	300x350
<input type="checkbox"/>	93363	Varilite Meridian™ PSV	350x350
<input type="checkbox"/>	93364	Varilite Meridian™ PSV	350x400
<input type="checkbox"/>	93411	Varilite Meridian™ PSV	350x450
<input type="checkbox"/>	93412	Varilite Meridian™ PSV	350x500
<input type="checkbox"/>	93414	Varilite Meridian™ PSV	380x380
<input type="checkbox"/>	93415	Varilite Meridian™ PSV	380x430
<input type="checkbox"/>	93413	Varilite Meridian™ PSV	400x350
<input type="checkbox"/>	93365	Varilite Meridian™ PSV	400x400
<input type="checkbox"/>	93366	Varilite Meridian™ PSV	400x450
<input type="checkbox"/>	93367	Varilite Meridian™ PSV	400x500
<input type="checkbox"/>	93368	Varilite Meridian™ PSV	430x430
<input type="checkbox"/>	93369	Varilite Meridian™ PSV	450x400
<input type="checkbox"/>	93370	Varilite Meridian™ PSV	450x450
<input type="checkbox"/>	93371	Varilite Meridian™ PSV	450x500
<input type="checkbox"/>	93571	Varilite Meridian™ PSV	450x550
<input type="checkbox"/>	93410	Varilite Meridian™ PSV	500x400
<input type="checkbox"/>	93372	Varilite Meridian™ PSV	500x450
<input type="checkbox"/>	93373	Varilite Meridian™ PSV	500x500
<input type="checkbox"/>	93416	Varilite Meridian™ PSV	550x450
<input type="checkbox"/>	93417	Varilite Meridian™ PSV	550x500
<input type="checkbox"/>	93418	Varilite Meridian™ PSV	600x450
<input type="checkbox"/>	93419	Varilite Meridian™ PSV	600x500



For wheelchair users with challenging, asymmetric skin protection, positioning and postural support requirements. Modifiable.

Varilite Pro-Form NX™  
Dual chamber

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)
<input type="checkbox"/>	93296	Varilite Pro-Form NX™, dual chamber	300x350
<input type="checkbox"/>	93297	Varilite Pro-Form NX™, dual chamber	350x350
<input type="checkbox"/>	93298	Varilite Pro-Form NX™, dual chamber	350x400
<input type="checkbox"/>	93113	Varilite Pro-Form NX™, dual chamber	380x380
<input type="checkbox"/>	93114	Varilite Pro-Form NX™, dual chamber	380x430
<input type="checkbox"/>	93299	Varilite Pro-Form NX™, dual chamber	400x400
<input type="checkbox"/>	93300	Varilite Pro-Form NX™, dual chamber	400x450
<input type="checkbox"/>	93301	Varilite Pro-Form NX™, dual chamber	400x500
<input type="checkbox"/>	93115	Varilite Pro-Form NX™, dual chamber	430x430
<input type="checkbox"/>	93302	Varilite Pro-Form NX™, dual chamber	450x400
<input type="checkbox"/>	93303	Varilite Pro-Form NX™, dual chamber	450x450
<input type="checkbox"/>	93304	Varilite Pro-Form NX™, dual chamber	450x500
<input type="checkbox"/>	93116	Varilite Pro-Form NX™, dual chamber	500x400
<input type="checkbox"/>	93117	Varilite Pro-Form NX™, dual chamber	500x450
<input type="checkbox"/>	93305	Varilite Pro-Form NX™, dual chamber	500x500



Non-adjustable skin protection cushion featuring VARILITE's unique Air-Release Device™.

Varilite Reflex™

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)
<input type="checkbox"/>	93491	Varilite Reflex™	300x300
<input type="checkbox"/>	93492	Varilite Reflex™	300x350
<input type="checkbox"/>	93493	Varilite Reflex™	350x350
<input type="checkbox"/>	93494	Varilite Reflex™	350x400
<input type="checkbox"/>	93495	Varilite Reflex™	350x450
<input type="checkbox"/>	93496	Varilite Reflex™	350x500
<input type="checkbox"/>	93497	Varilite Reflex™	380x380
<input type="checkbox"/>	93498	Varilite Reflex™	380x430
<input type="checkbox"/>	93499	Varilite Reflex™	400x350
<input type="checkbox"/>	93500	Varilite Reflex™	400x400
<input type="checkbox"/>	93501	Varilite Reflex™	400x450
<input type="checkbox"/>	93502	Varilite Reflex™	400x500
<input type="checkbox"/>	93503	Varilite Reflex™	430x430
<input type="checkbox"/>	93504	Varilite Reflex™	450x400
<input type="checkbox"/>	93505	Varilite Reflex™	450x450
<input type="checkbox"/>	93506	Varilite Reflex™	450x500
<input type="checkbox"/>	93507	Varilite Reflex™	500x400
<input type="checkbox"/>	93508	Varilite Reflex™	500x450
<input type="checkbox"/>	93509	Varilite Reflex™	500x500



Varilite Stratus

The Varilite Stratus Cushion is a great choice for people who have a moderate to low risk of tissue breakdown and moderate positioning needs. The cushion has a unique self-inflating valve, which allows air in and out of the cushion so it can contour to the user's shape.

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)
<input type="checkbox"/>	93047	Varilite Stratus	300x300
<input type="checkbox"/>	93048	Varilite Stratus	300x350
<input type="checkbox"/>	93049	Varilite Stratus	350x350
<input type="checkbox"/>	93050	Varilite Stratus	350x400
<input type="checkbox"/>	93110	Varilite Stratus	380x380
<input type="checkbox"/>	93111	Varilite Stratus	380x430
<input type="checkbox"/>	93051	Varilite Stratus	400x400
<input type="checkbox"/>	93052	Varilite Stratus	400x450
<input type="checkbox"/>	93053	Varilite Stratus	430x430
<input type="checkbox"/>	93054	Varilite Stratus	450x400
<input type="checkbox"/>	93055	Varilite Stratus	450x450
<input type="checkbox"/>	93112	Varilite Stratus	450x500
<input type="checkbox"/>	93056	Varilite Stratus	500x450
<input type="checkbox"/>	93057	Varilite Stratus	500x500



Varilite Zoid™  
PSV

Low-profile, adjustable skin protection performance cushion for active wheelchair users

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)
<input type="checkbox"/>	93031	Varilite Zoid™ PSV	300x300
<input type="checkbox"/>	93032	Varilite Zoid™ PSV	300x350
<input type="checkbox"/>	93033	Varilite Zoid™ PSV	350x350
<input type="checkbox"/>	93034	Varilite Zoid™ PSV	350x400
<input type="checkbox"/>	93100	Varilite Zoid™ PSV	350x450
<input type="checkbox"/>	93101	Varilite Zoid™ PSV	380x380
<input type="checkbox"/>	93102	Varilite Zoid™ PSV	380x430
<input type="checkbox"/>	93035	Varilite Zoid™ PSV	400x400
<input type="checkbox"/>	93036	Varilite Zoid™ PSV	400x450
<input type="checkbox"/>	93103	Varilite Zoid™ PSV	400x500
<input type="checkbox"/>	93104	Varilite Zoid™ PSV	430x430
<input type="checkbox"/>	93105	Varilite Zoid™ PSV	450x400
<input type="checkbox"/>	93037	Varilite Zoid™ PSV	450x450
<input type="checkbox"/>	93106	Varilite Zoid™ PSV	450x500
<input type="checkbox"/>	93107	Varilite Zoid™ PSV	500x400
<input type="checkbox"/>	93108	Varilite Zoid™ PSV	500x450
<input type="checkbox"/>	93109	Varilite Zoid™ PSV	500x500



Varilite Junior

Ideal for children with a high degree of tissue breakdown and symmetric positioning needs

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)
<input type="checkbox"/>	93588	Varilite Junior	250x250
<input type="checkbox"/>	93599	Varilite Junior	275x325
<input type="checkbox"/>	93600	Varilite Junior	300x300

## 5 Backrest

### Backrest bracket

CODE	DESCRIPTION	SEAT WIDTH (mm)
<input type="checkbox"/> 80298	Backrest bracket adjustable	seat width adjustable in steps 380 - 400 - 450 - 500 - 550 mm
<input type="checkbox"/> 80298-1	Backrest bracket adjustable JUNIOR	seat width adjustable in steps 250 - 300 - 350 mm
<input type="checkbox"/> 80355	Backrest bracket adjustable Extra long tubes for mounting other back systems	seat width adjustable in steps 380 - 400 - 450 - 500 - 550 mm



80298



80298-1  
Junior

### Backrest system



Icon™ Back System  
Low

Combines innovative VariLock™ hardware with a sleek, low profile look



Icon™ Back System  
Mid

Provides support and positioning to improve trunk control and decrease fatigue



Icon™ Back System  
Tall

Well suited for users needing additional support or with a tilt-recline system



Icon™ Back System  
Deep

Ideal for users needing maximum support

CODE	DESCRIPTION	SEAT WIDTH (mm)
<input type="checkbox"/> 93560	Icon™ Back System Low	350
<input type="checkbox"/> 93561	Icon™ Back System Low	380
<input type="checkbox"/> 93562	Icon™ Back System Low	400
<input type="checkbox"/> 93563	Icon™ Back System Low	430
<input type="checkbox"/> 93564	Icon™ Back System Low	450
<input type="checkbox"/> 93565	Icon™ Back System Low	500
<input type="checkbox"/> 93587	Icon™ Back System Mid	250
<input type="checkbox"/> 93544	Icon™ Back System Mid	300
<input type="checkbox"/> 93545	Icon™ Back System Mid	350
<input type="checkbox"/> 93546	Icon™ Back System Mid	380
<input type="checkbox"/> 93547	Icon™ Back System Mid	400
<input type="checkbox"/> 93548	Icon™ Back System Mid	430
<input type="checkbox"/> 93549	Icon™ Back System Mid	450
<input type="checkbox"/> 93550	Icon™ Back System Mid	500
<input type="checkbox"/> 93551	Icon™ Back System Mid	550
<input type="checkbox"/> 93552	Icon™ Back System Mid	600
<input type="checkbox"/> 93553	Icon™ Back System Tall	300
<input type="checkbox"/> 93554	Icon™ Back System Tall	350
<input type="checkbox"/> 93555	Icon™ Back System Tall	380
<input type="checkbox"/> 93556	Icon™ Back System Tall	400
<input type="checkbox"/> 93557	Icon™ Back System Tall	430
<input type="checkbox"/> 93558	Icon™ Back System Tall	450
<input type="checkbox"/> 93559	Icon™ Back System Tall	500
<input type="checkbox"/> 93586	Icon™ Back System Deep	250
<input type="checkbox"/> 93538	Icon™ Back System Deep	350
<input type="checkbox"/> 93539	Icon™ Back System Deep	380
<input type="checkbox"/> 93540	Icon™ Back System Deep	400
<input type="checkbox"/> 93541	Icon™ Back System Deep	430
<input type="checkbox"/> 93542	Icon™ Back System Deep	450



---

<input type="checkbox"/> 93543	Icon™ Back System Deep	500
--------------------------------	------------------------	-----

---



### Backrest adjustments and options

CODE	DESCRIPTION	
<input type="checkbox"/> 80363	Mechanical backrest adjustment, quick release	
<input type="checkbox"/> 80577	Mechanical backrest adjustment, gas spring	
<input type="checkbox"/> 97235	Electric backrest	95°-135°
<input type="checkbox"/> 81137	Electric backrest with extender bracket actuator	95°-150°



80363



80577

### Lateral support

CODE		DESCRIPTION					
for seat width							
30cm	35-38cm	40-43cm	45cm	50-60cm			
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Varilite lateral support fixed, set	Choose underneath Cushion for thigh/lateral support for this.	
93254	93255	93256	93257	93258			
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Varilite lateral support swing away, set	short Choose underneath Cushion for thigh/lateral support for this.	
93239	93240	93241	93242	93243			
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Varilite lateral support swing away, set	medium Choose underneath Cushion for thigh/lateral support for this.	
93249	93250	93251	93252	93253			
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Varilite lateral support swing away, set	long Choose underneath Cushion for thigh/lateral support for this.	
93244	93245	93246	93247	93248			
<input type="checkbox"/>	93262	Cushion for thigh/lateral support, pcs 7,5x10 cm, pcs.					
<input type="checkbox"/>	93263	Cushion for thigh/lateral support, pcs 10x13 cm, pcs.					
<input type="checkbox"/>	93264	Cushion for thigh/lateral support, pcs 13x15 cm, pcs.					



Fixed



Swing away short



Swing away medium



Swing away long

---

### Head support

---

	CODE	DESCRIPTION
<input type="checkbox"/>	80337	Head support (225 x 120 mm) incl. bracket (80288)
<input type="checkbox"/>	80829	Head support medium (175 x 95 mm) incl. bracket (80288-2)
<input type="checkbox"/>	80267	Head support with support pad

---



80267

---

### Thigh support

---

	CODE	DESCRIPTION
<input type="checkbox"/>	93262 80328	Thigh support, incl. bracket, cushion 75 x 100 mm, set
<input type="checkbox"/>	93263 80328	Thigh support, incl. bracket, cushion 100 x 130 mm, set
<input type="checkbox"/>	93264 80328	Thigh support, incl. bracket, cushion 130 x 150 mm, set

---



## Belts & support

CODE					DESCRIPTION
<input type="checkbox"/>	93231				Varilite chest belt
	<b>XS</b>	<b>S</b>	<b>M</b>	<b>L</b>	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Varilite 2-point hip belt
	93221	93222	93223	93224	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Varilite 4-point hip belt
	93225	93226	93227	93228	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Varilite Band-It, Ankle, set
	93235	93236	93237	93238	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Varilite contoured chest harness
	93272	93232	93233	93234	
<input type="checkbox"/>	80765				Bracket for chest harness
<input type="checkbox"/>	97157				Belt, automatic roll-up



Chest belt



2-point hip belt



4-point hip belt



Ankle



Harness

## 6 Arm support

CODE	DESCRIPTION	
<input type="checkbox"/>	80289-2 Arm support flip-up with front support lockable, set	height adjustable from 185 to 285 mm
<input type="checkbox"/>	80302 Arm support flip-up, set	height adjustable from 185 to 285 mm
<input type="checkbox"/>	80743 Stop for Flip up arm support, pcs	prevents arm support from going further back
<input type="checkbox"/>	80347 Arm support bracket rigid, set	adjustable from 100 to 280 mm from the seat plate

## Options for armsupport

CODE	DESCRIPTION			
<input type="checkbox"/> 80320	Arm pad cushion, right, short 320mm		length 320mm	
<input type="checkbox"/> 80335	Arm pad cushion, right, long 430mm		length 430mm	
<input type="checkbox"/> 80336	Arm pad cushion, left, short 320mm		length 320mm	
<input type="checkbox"/> 80321	Arm pad cushion, left, long 430mm		length 430mm	
<input type="checkbox"/> 80651	Aluminium edge for arm support cushion short			
<input type="checkbox"/> 80652	Aluminium edge for arm support cushion long			
<input type="checkbox"/> 80371	Arm support bracket, set	<input type="checkbox"/> VL pad 75x100mm	<input type="checkbox"/> VL pad 100x130mm	<input type="checkbox"/> VL pad 130x150mm
<input type="checkbox"/> 80830	Swing bracket for arm pad right		for armrest flip-up 80302	
<input type="checkbox"/> 80830-1	Swing bracket for armpad left		for armrest flip-up 80302	
<input type="checkbox"/> 80994	Purse -leather bag for arm support			
<input type="checkbox"/> 20248	Purse -pocket for arm support			



80320/80336



80335/80321



80651/80652



80830/80830-1



80371



80994



20248

## 7 Leg support and Footplate

### Two-piece leg support

CODE	DESCRIPTION		
<input type="checkbox"/> 80700-1	Manual elevating leg support, right		
<input type="checkbox"/> 80701-1	Manual elevating leg support, left		
<input type="checkbox"/> 80628	Electric elevating leg support, right		
<input type="checkbox"/> 80629	Electric elevating leg support, left		
<input type="checkbox"/> 80716	Amputation elevating leg support right		
<input type="checkbox"/> 80717	Amputation elevating leg support left		
<input type="checkbox"/> 80721	Knee pad standard right		for electric separate leg support, right 80628
<input type="checkbox"/> 80730	Knee pad standard left		for electric separate leg support, left 80629
<input type="checkbox"/> 80631-3	Lower tube for leg support straight + 100mm		
<input type="checkbox"/> 80631-2	Lower tube for leg support curved + 100mm		



80700-1/80701-1



80628/80629



80716/80717



80721 / 80730

### Central leg support with single footplate

CODE	DESCRIPTION	
<input type="checkbox"/> 95083	Center manuel leg support adjustable	
<input type="checkbox"/> 80357	Center power leg support with calf pads	mechanically lenght adjustment (via screws) 340-430mm + 150mm electrical length adjustment (via controller)
<input type="checkbox"/> 81124	Center power leg support with calf pads for sepearate footplates	mechanically lenght adjustment (via screws) 340-430mm + 150mm electrical length adjustment (via controller)
<input type="checkbox"/> 80571	Center power leg support short with calfpads	mechanically lenght adjustment (via screws) 270-310mm + 80mm electrical length adjustment (via controller)
<input type="checkbox"/> 81126	Center power leg support short with calfpads for separate footplates	mechanically lenght adjustment (via screws) 250-310mm + 80mm electrical length adjustment (via controller)
<input type="checkbox"/> 81130 81017	Center power leg support short with calfpads incl. electric flip-up footplate	mechanically lenght adjustment (via screws) 270-390mm + 80mm electrical length adjustment (via controller)
<input type="checkbox"/> 81128 81017	Center power leg support with calf pads incl. electric flip-up footplate	mechanically lenght adjustment (via screws) 340-430mm + 150mm electrical length adjustment (via controller)



95083



81017  
Footplate flip-up  
electric



80357

### Footplate for central leg support

CODE	DESCRIPTION
<input type="checkbox"/> 80373-2	Footplate moulded standard complete with edge
<input type="checkbox"/> 80539-1	Footplate with cutout moulded complete with edge
<input type="checkbox"/> 80818	Footplates divided



80373-2



80539-1



80818

### Options for central leg support

CODE	DESCRIPTION
<input type="checkbox"/> 80690	Leather pad for center leg support, knee protection



80690

## 8 Lighting

CODE	DESCRIPTION
<input type="checkbox"/> 97223	Active LED lighting
<input type="checkbox"/> ---	Without lighting

## 9 Control module and accessories

CODE	DESCRIPTION		
<input type="checkbox"/> 95174	Controller, JSM, LED Entry level	w/o light	Power Module 120A
<input type="checkbox"/> 95253	Controller, JSM-L LED	for light	Power Module 120A
<input type="checkbox"/> 95172-1	Controller, CJSM, SW (direct actuator)	w/o light	Power Module 120A
<input type="checkbox"/> 95252-1	Controller, CJSM-L-SW(direct actuator)	for light	Power Module 120A
<input type="checkbox"/> 95501	Controller CJSM2-L	infrared function	Power Module 120A
<input type="checkbox"/> 95524	Controller, CJSM2, Bluetooth	infrared & bluetooth funct.	Power Module 120A



95174/95253

JSM LED



95172-1/95252-1

CJSM



95501/95524

CJSM-2

## Control module bracket

CODE	DESCRIPTION
<input type="checkbox"/> 80407	Right swingaway bracket for control module right
<input type="checkbox"/> 80408	Left swingaway bracket for control module left
<input type="checkbox"/> 80838	Drop-Down bracket universal



80407/80408



80838

### Control module knob

	CODE	DESCRIPTION
<input checked="" type="checkbox"/>	95153	Standard
<input type="checkbox"/>	95154	Mushroom M48
<input type="checkbox"/>	95155	Ball 38 mm
<input type="checkbox"/>	95299	Ball 32mm
<input type="checkbox"/>	95306	Chin C38
<input type="checkbox"/>	95298	Cone S30
<input type="checkbox"/>	95305	Stick S80
<input type="checkbox"/>	95301	T-bar
<input type="checkbox"/>	95300	Softball 70mm
<input type="checkbox"/>	Select color	Golfball



95153



95154



95155/95299



95306



95298



95305



95301



95300



White  
 80653



Orange  
 80653-4



Pink  
 80653-2



Yellow  
 80653-3



Football  
 80653-7



Football  
 80653-8

### Special operating modules

	CODE	DESCRIPTION	
<input type="checkbox"/>	97342	TA hand control for CJSM, incl. bracket (80678)	
<input type="checkbox"/>	97343	TA hand control for CJSM2, incl. bracket (80678)	
<input type="checkbox"/>	95213	Attendant module incl. connector block and canbus	
<input type="checkbox"/>	80903	1 hand Attendant bracket swing	without controller 95213
<input type="checkbox"/>	80517	1 hand attendant bracket IW / iQ	without controller 95213
<input type="checkbox"/>	80533	TA 2-hand attendant control, with emergency stop	without controller 95213



97342/97343

## 10 Accessories

	CODE	DESCRIPTION	
<input type="checkbox"/>	80163	Tray, adult, incl. bracket (80518)	
<input type="checkbox"/>	80164	Tray, junior, incl. bracket (80518)	
<input type="checkbox"/>		Half tray, on request and customized	
<input type="checkbox"/>	80175	Bag hooks	
<input checked="" type="checkbox"/>	97285	TA flight kit	Flight kit is designed to disconnect the batteries easily. This is an requirement for most airlines



80163/80164



80175