



THIS IS THE BENCHMARK:

Double-Drum test & crash test
200,000 load changes without overload

THIS IS HOW WE TEST:

Double-Drum test & crash test
and additional MEYRA test procedures
200,000 - 1,000,000 LOAD CHANGES
with 10 % overload + additional stress test

THEREFORE MEYRA®!

WHY MEYRA®?

GENUINE GERMAN QUALITY

MORE THAN
20,000
PRODUCTION UNITS PER YEAR
20,000 m²
COVERED PRODUCTION AREA

OVER
50,000 m²
SITE AREA

... AND MOREOVER, A MULTIPLE
AWARD-WINNING HEALTH CARE SYSTEM
FOR EMPLOYEES

2 RESILIENT HIGHEST TEST REQUIREMENTS

TEST TRACKS & CRASH TESTS

All wheelchairs are extensively tested on a test track. A wheelchair is only released for delivery after successful completion of the internal quality audits. In addition, crash tests are regularly performed on new models to check the resilience of our wheelchairs. Findings from possible deformations are used to optimize a design permanently.

MEYRA ENDURANCE TEST

MEYRA tests beyond the standard of legal standards. Our quality management system is certified in accordance with DIN EN ISO 13485:2012. The requirements of the standards are tested annually and confirmed by a test report. In addition, we carry out stress and reliability tests with increased stress beyond the statutory standard. This is unique and ensures adherence to our high quality demands on our wheelchairs.

MEYRA GmbH

Meyra-Ring 2
D-32689 Kalletal-Kalldorf

Tel +49 5733 922 - 0
Fax +49 5733 922 - 9311

info@meyra.de

www.meyra.de



FOUR GOOD REASONS





5 YEARS GUARANTEE*

Extremely stable - extremely durable!
Now available: Take advantage of the
extended warranty for the entire iCHAIR
family and Optimus.



1

DURABLE OUR PROMISE OF HIGH QUALITY

WELDING TECHNOLOGY REQUIRES...

a lot of effort, which is worth it! For the designing of our framework, we rely solely on the so-called inert gas welding processes MAG and TIG. This lets us to provide an optimised surface. For guaranteed durability and low maintenance requirements.

ESPECIALLY HIGH-QUALITY SURFACES

With modern electroplating technique, we lay the foundation for the corrosion resistance and guarantee an exceptionally high durability of our frame. Electrostatic powder coating ensures high scratch resistance. We attach particular importance to high environmental compatibility. The use of single-origin materials, clear material identification as well as the monitoring at each and every stage of the process guarantee compliance with our environmental standards.

* 5-year guarantee on the frame. See MEYRA guarantee conditions.

3

SOLUTION-ORIENTED ACCORDING TO DISEASE PATTERNS

ALWAYS WELL ADVISED

Different disease patterns mean individual requirements for the functionality of a wheelchair - a demanding task for the specialist trade in consulting. We have set ourselves the target of providing you with optimal support. Whether paediatrics, palliative medicine, geriatrics, neurology, orthopedics / traumatology or internal medicine - we want to ensure the best possible care

with our lucid assignment table. You may find a list of our products in our product catalog, with which disease patterns are certain wheelchairs used and in which medical specialty are they used. For good advice and satisfied customers.

**MORE INFORMATION
IN THE PRODUCT CATALOG**



4

MODULAR EASY SPAREPARTS SUPPLY

High efficiency through maximum parts compatibility - this idea lies behind the MEYRA modular system. This is already an integral part of the EUROCHAIR and iCHAIR families. The modular system has also largely established itself in other product areas. The equality of the single side parts installed in our wheelchairs also guarantees a highest level of flexibility and excellent user-friendliness.

**MORE INFORMATION
IN K.I.S.S.S., PRODUCT CATALOGUE OR VIDEO**



> Direct to video