



All seat dimensions are always **measured** without seat cushion!

Body measurements

Please refer to the sketch in the enclosed information!

MEASUREMENT	ON USER	ON WHEELCHAIR
01	Seat depth in mm From the hollow of the knee to the back of the buttocks	From the back to the front edge of the seat pan
02	Seat width in mm Widest part of the pelvis when sitting	Width of the seat pan or distance between the side panels
03	Chest width in mm Distance between the armpits	
04	Chest depth in mm Under the armpit, from the back cushion to the front of the thorax/chest	
05	Elbow height in mm Distance between the upper edge of the seat-pan and the lower edge of the elbow (90° elbow angle)	
06	Armpit height in mm Distance between the upper edge of the seat pan and the armpits	
07	Back height in mm Seat pan to upper edge of shoulder	Upper edge of seat pan to upper edge of back strap
08	Head height in mm Distance between the upper edge of the seat pan and the highest point of the head	
09	Forearm length in mm Distance between the elbow and the bent wrist	
10	Lower leg length in mm From the hollow of the knee to the lower edge of the heel (incl. shoe)	Upper edge of footplate to upper edge of seat pan