

TECHNICAL DATA

Dimensions in mm, weights in kg. Wheels in mm, speeds in km / h, tolerances + - 10mm, + - 2°, data subject to design changes

CL 409 1.254

Length without basket in mm	1.200
Width in mm	610
Seat height in mm, adjustable	420 – 480
Seat depth in mm	410
Back height in mm	450
Seat width in mm, adjustable	450 – 670
Transport height (with steering column folded and without seat) in mm	600
Number of wheels	4
Diameter steering wheel in mm	260
Diameter drive wheel in mm	260
Ground clearance in mm	100
Obstacle overcoming, max in mm	60
Turning circle (radius) in mm	1.500
Motor continuous power in W	210 W
Motor peak power in W	850 W
Batteries in Ah	28,1 Ah (5 h) / 35 Ah (20 h)
Charger in A	4 A
Electronics in A	70 A
Max. range (under test conditions) in km	30
Percentage allowed in ascents / descents	9° (16 %)
Empty weight in kg with batteries	91,5
Permissible user weight in kg	136
Permitted total weight in kg	230
Max. speed in km / h	6

MEYRA GmbH

📍 Meyra-Ring 2
D-32689 Kalletal-Kalldorf

☎ Tel +49 5733 922 - 555
Fax +49 5733 922 - 9311

✉ info@meyra.de

www.meyra.com

by www.bwc.gmbh



ORTOPEDIA
EINE MARKE VON MEYRA

CL 409
The economic
electromobile

E-POWER

A DIFFERENT WAY FOR A CITY STROLL



The new CL 409 is a reliable companion for everyday life – whether shopping in the shopping center, sightseeing in the city or a trip to the park. Its seating system with infinitely adjustable back angle and high seating comfort is unique in its product class. It convinces its users with a comfortable seating position even on longer distances. The dashboard is clearly structured and ensures simple and intuitive handling – so nothing stands in the way of an endless city tour!

Latest economical LED technology



Simple dashboard



Sprung comfort seat system, large legroom and stepless seat adjustment



Extremely robust platform

Modern design with pleasing shapes in trendy colors



Aluminum rims