CL510

1.264 Stock version

Allowed user weight

165 kg (including payload)



ORTOPEDIA

EINE MARKE VON MEYRA



LED technology



1.600 Watt motor



Extremely robust cassette-type frame



Spring-loaded comfort seating system

The powerful reliable

ADAPTATION

- The angle of the backrest, back extension and height of the seat and armrests are adjustable
- The seat width can be stepless adjusted via the armrests
- The seat has adjustment rails and can be rotated through 360°

OPERATION

- Very easy to use
- Manual smooth steering
- Driving lever for thumb operation and finger operation
- Mechanical release of the spring-loaded brake
- Self-explanatory control panel with speed preselection options

● ORDER ❍ COST ESTIMATE ❍ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

ORTOPEDIA
EINE MARKE VON MEYRA



Standard configuration

- alloy wheels
- 2 x rearview mirrors
- Reversible and lockable handlebar
- spring-loaded comfort seat
- Height- und angle-adjustable backrest, rotatable about 360°
- seat adjustment
- Swing-up armrests
- Double lever for thumb operation
- Lighting including turn signal, hazard warning light and brake light
- self-explanatory Display
- Battery charger
- horn
- Speed Preselection

Allowed user weight

165 kg (including payload)

Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply

Valid from 01.01.2026 – 31.12.2026.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

SCOOTER

Technical specifications

Length in mm	1200
Width in mm	610
Seat height in mm	420 - 480
Seat depth in mm	410
Seat width adjustable in mm	450 - 670
Back support height in mm	450
Transport height (with folded-down steering column and without seat)	600
Castor wheels diameter	260
Drive wheels diameter	260
Ground clearance in mm	100
Max. obstacle height in mm	60
Turning radius in mm	1500
Number of wheels	4
Suspension seat	
Continuous motor performance	400 W
Peak motor performance	1600 W
Batteries	2x 12V 50AH
Battery charger	6A
Electronic	120A
Range	30 km
Permitted inclination/slopes	9° / 16%
Empty weight from	93 kg
Permissible user weight in kg	165 kg
Permitted total weight in kg	260 kg

Please enter your quantities in the boxes

01 City liner

Item	DESCRIPTION
9126406	6 km/h- Dolphin grey metallic
9126416	6 km/h- Dolphin grey metallic - Handbrake system for Bus transport
9126410	10 km/h- Dolphin grey metallic
9126420	10 km/h- Dolphin grey metallic - Handbrake system for Bus transport

20 Options

Item	DESCRIPTION
3086642	Crutch holder
3075902	Body sack
1058996	Shopping basket additionally, not in combination with walking aid holder Code 3086642
1058988	Rear marking plate
1061128	TÜV certificate
12642001	*Sitzbreite Plus Variante 450 - 700 mm HMV-Nr.: 18.51.05.2001