TA IQ MWD

2.518

Allowed user weight

140 kg Allowed user weight



MySeat with comfort and flexibility

CHARACTERISTICS

- New seat adjustment module (tilt module) ensures
  greater stability and comfort for the user
- Biomechanical solutions for the backrest, armrests and leg rests with the MYSeat system
- The new design improves product quality and reliability, as well as service times for specialist retailers
- The best of the MEYRA Group in a single product:
  - TA iQ Powerbase
  - MYSeat System
  - NETTI Seating & Positioning
  - VARILITE Seating & Positioning

FIELD OF USE

- Extremely small turning circle for effortless manoeuvring in the tightest of spaces
- Outstanding agility indoors and outdoors for maximum flexibility in everyday life
- Six wheels ensure constant contact with the ground
  for excellent stability and a smooth ride
- Impressive off-road capability for greater independence
  across a range of terrain conditions

THERAPEUTIC BENEFIT

- Improving functional independence in everyday life and expanding the range of movement to enable independent use of shelves, work surfaces and control panels
- Promoting social interaction on an equal footing to
  improve communication and participation
- Improving the ability to carry out activities of daily living (ADL) independently, and thereby social participation
- TA IQ is suitable for
  - Muscular dystrophy
  - Multiple sclerosis (MS)
  - Cerebral palsy
  - Spinal cord injuries
  - Amyotrophic lateral sclerosis (ALS)
  - Traumatic brain injuries
  - Parkinson's disease
  - Severe rheumatoid arthritis

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

Order quantity:

MEYRA



Price basic model ***16.990,00 €***

Allowed user weight
140 kg Allowed user weight

Hints
Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01.01.2026 – 31.12.2026
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

USER DATA

First name

Name

Street

ZIP CODE

City

Country

Birthday

Diagnosis 1

Diagnosis 2

Tel.

Mobil

E-mail

Order number

Health insurance fund / cost unit

Specialist trade Contact person/phone.

Supply specialist MEYRA

Approval info from week

ANATOMICAL DATA (IN MM)

1) Seat depth

2) Seat width

3) Chest width

4) Chest depth

5) Elbow height

6) Armpit height

7) Back height


8) Head height

9) Forearm length

10) LLL

Body height

Weight (kg)



Remarks

❑ I accept the Terms and Conditions of Meyra Distribution GmbH (valid as of January 2023)
https://www.meyra.com/gtc/

Date

City

Signature

TAS

Technical specifications

Overall length (measured at 0° seat angle) in mm	1190 - 1350
Total width in mm	600 - 770
Overall height in mm	910 - 1180
Permissible user weight in kg	140
Permitted total weight in kg	320
Transport empty weight in kg	100 kg without battery / 160 kg with 80Ah battery
Motor power	2 x 350 Watt
Speeds in km/h	6 / 10 / 12
Range in km	35 - 40
Obstacle height in mm	max. 80
Permitted inclination/slopes	10° (18 %)
Electrically adjustable backrest angle	80° - 140° / 90° - 150° / 90° - 170° (Code 25 / Code 2947 / Code 2948)
Mechanical backrest angle	80° - 120° (Code 401)
Seat stroke in mm	300
Seat camber electric	0° - 45°
Minimum turning radius in mm	700
Ground clearance in mm	60 (12,5" drive wheel), 70 (14" drive wheel)

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

01 Chassis / engine / frame options

Chassis

	CODE	DESCRIPTION	Drive type	EUR €
❑	2528	IQ MWD	Mid-wheel drive	0,00

Engine

	CODE	DESCRIPTION	RANGE	MOTOR POWER	MAX. USER WEIGHT	EUR €
❑	2115	6 km/h	max. 40 km	2 x 350W 4-POL	140 kg	0,00
❑	2117	10 km/h	max. 35 km	2 x 350W 4-POL	140 kg	490,00
❑	2120	12 km/h	max. 35 km	2 x 350W 4-POL	140 kg	690,00

ATTENTION:
For 10 km/h and 12 km/h motorised versions:
in Germany, only in conjunction with TÜV certificate CODE 4772
plus rear marking film CODE 2676
or rear marking plate CODE 676

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

02 Wheels and tyres

Drive wheels

		FSH (mm)			
		(C-2511)	(C-2512)	(C-2511)	(C-2512)
Code	**Ø Rim**	**385**	**405**	**395**	**415**
815	12"	❑	❑		
460	14"			❑	❑

⚠ *12" drive wheel CODE 815 not compatible with Dahl docking station adapter set CODE 2721*
VSH 385 mm / 395 mm CODE 2511 compatible only with:
armrest CODE 24
central leg support CODE 5100
central footplate CODE 5109

Drive wheels Tyres

	CODE	DESCRIPTION	Tyres	EUR €
❑	482	Tread pattern	Pneumatic tyres	0,00
❑	905	Tread pattern	Puncture-proof tyres	228,00

Steering wheels / Tyres

		Tyres	
		(C-146)	(C-900)
Code	**Ø Lightwheel**	**Pneumatic tyres**	**Puncture-proof tyres**
157	8" 200 x 50 mm	❑	❑ 228,00
163	8,5" 220 x 65 mm	❑	❑ 228,00

⚠ *8,5" steering wheel CODE 163 not compatible with 12' drive wheel CODE 815*

TAS

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Rim colour

	CODE	DESCRIPTION	EUR €
❑	200	Black	0,00
❑	201	Silver-black ⚠ *Not compatible with 8.5" steering wheel CODE 163* *Not compatible with 12" drive wheel CODE 815*	0,00

Tyre colour

	CODE	DESCRIPTION	EUR €
❑	909	Grey	0,00
❑	908	Black ⚠ *Not compatible with 12" drive wheel CODE 815*	0,00

Mud guard

	CODE	DESCRIPTION	EUR €
❑	816	Mud guard for Castor wheels ⚠ *Not compatible with 8.5" steering wheel CODE 163*	146,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

03 Seat

Seat width/depth

		Seat depth (mm)					
		(C-615)	(C-788)	(C-619)	(C-4105)	(C-791)	(C-955)
Code	**Seat width**	**400**	**430**	**460**	**490**	**520**	**550**
38	380 mm	❑	❑	❑	❑	❑	❑
43	430 mm	❑	❑	❑	❑	❑	❑
48	480 mm	❑	❑	❑	❑	❑	❑
53	530 mm	❑	❑	❑	❑	❑	❑

⚠ *Seat depth 550 mm CODE 955 not compatible with electric seat tilt CODE 4118*

Seat Cushions

	CODE	DESCRIPTION	EUR €
❑		without seat cushion	0,00
❑	948	iCHAIR MEYLIFE contour seat cushion	412,00
❑	7206	Netti Uno	240,00
❑	7203	Netti Sit cushion	460,00

⚠ *NETTI cushion (not included) with a seat depth of 550 mm CODE 955*

TAS

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

VARILITE seat cushion

Please use the separate order form 'VARILITE Seating System' when selecting VARILITE seat cushions.
The 'VARILITE Seating System' order form can be found on the MEYRA website under model 2.518 TA IQ.
When placing an order, please enclose the completed VARILITE order form with this order form.

⚠ *Only in conjunction with the round-tube backrest CODE 608*

Cover for MYSeat contour cushion

	CODE	DESCRIPTION	EUR €
❑	237	Textile cover black	0,00
❑	254	Synthetic leather	434,00

⚠ *is only used on contour cushions 948/949*

Seat tilt / seat lift

	CODE	DESCRIPTION
❑	4118	Electric tilt, adjustable from 0° to 45°
❑	27	Seat lift 300 mm

⚠ *Both options are included in the basic model*

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

04 Back systems

Round tube back, back cushion and height

		Back height (mm)	
		(C-921)	(C-922)
Code	**Description**	**520**	**570**
608	Round tubes, 22 mm, prepared for external backrest systems	❑	❑
608/735	Strap back without polster	❑	❑
608/736	Strap back with polster, black	❑	❑
608/7201	Netti Uno Back, 3D, breathable	❑ 396,00	❑ 396,00
608/7209	Netti Smart, 3D, breathable	❑ 394,00	❑ 394,00
608/7207	Netti Stabil, 3D, breathable	❑ 648,00	❑ 648,00
608/7210	Netti Super Stabil, 3D, breathable	❑ 744,00	❑ 744,00
608/7205	Netti Kyphotic, 3D, breathable	❑ 882,00	❑ 882,00

VARILITE Back cushion

To select VARILITE backrest pads, please use the separate order form 'VARILITE Seating System'.
The 'VARILITE Seating System' order form can be found on the MEYRA website under model 2.518 TA IQ.
When placing your order, please enclose the completed VARILITE order form with this order form.

⚠ *Only in conjunction with the round-tube backrest CODE 608*

Back adjustments for round tube backs

	CODE	DESCRIPTION	EUR €
❑	401	Mechanical back recline for round tube backs	0,00
❑	25	Electric back recline for round tube backs	1.650,00

☐ Please select by ticking the relevant boxes (S) Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

Central back/contour back

		Back height (mm)			
		(C-921)	(C-922)	(C-923)	(C-928)
Code	**Description**	**520**	**570**	**620**	**670**
949/2038/2859	Contour back 380mm wide	❑ 372,00	❑ 372,00	❑ 372,00	❑ 372,00
949/2043/2859	Contour back 430mm wide	❑ 372,00	❑ 372,00	❑ 372,00	❑ 372,00
949/2048/2859	Contour back 480mm wide	❑ 372,00	❑ 372,00	❑ 372,00	❑ 372,00
949/2053/2859	Contour back 530mm wide	❑ 372,00	❑ 372,00	❑ 372,00	❑ 372,00

⚠ *CODE 921 and CODE 922 not compatible with electric biomechanical backrest adjustment for central backrest CODE 2948*
CODE 928 not compatible with electric backrest adjustment for central backrest CODE 2947

Back adjustments for central back

	CODE	DESCRIPTION	Degree angle		EUR €
❑	2946/2947	Electric backrest recline for central back	150°	infinitely variable	2.590,00
		⚠ *Not compatible with biomechanical armrest CODE 324*			
❑	2946/2948	Electric biomechanical back adjustment for central back	170°	infinitely variable	3.158,00
		⚠ *Only in conjunction with armrest swivel-up, angle-adjustable, biomechanical CODE 324*			

Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Additional options/positioning

	CODE	DESCRIPTION	POSITIONING	EUR €
❑	241	Multifunctional headrest with an easy-care cover		490,00
❑	242	Multifunctional headrest with 3D mesh cover		490,00
❑	937	Headrest extension		84,00
❑	4275	Mounting rails, accessories on the back *Not compatible with round tube backrest CODE 608*	Contour	332,00
❑	2957	Swing-away chest laterals *Only in conjunction with mounting rails CODE 4275*	Contour	590,00
❑	797	Side bolster for seat cushion (Seat Cushion side support)	Seat pan	77,00
❑	4933	Thigh Guides continuously adjustable in width and length	Seat pan	620,00
❑	814	Push handle continuous for central back *Not compatible with round tube backrest CODE 608*	Contour	298,00
❑	2958	Depth extension for round tube backs	Round tube	210,00
❑	2311	Back tubes outrigged 25 mm per side for round tube back Code 608	Round tube	316,00
❑	2313	Back tubes Inrigged 25 mm per side for round tube back Code 608	Round tube	316,00

05 Armrests

1 = removable 2 = height-adjustable 3 = width-adjustable 4 = depth-adjustable 5 = height-adjustable 6 = biomechanical

	CODE	DESCRIPTION	FUNCTIONS	EUR €
❑	106	Single post Armrest removable, height and depth adjustable *Not in combination with the electric biomechanical backrest adjustment CODE 2948* *Not in combination with seat height A CODE 2511*	1,2,3,4	0,00
❑	24	Armrest Flip up, angle-adjustable *Not i.c. with elect. biomechanical backrest adjustment CODE 2948*	2,3,4,5	640,00
❑	324	Armrest Flip up, angle-adjustable, biomechanical *Not compatible with electric backrest adjustment CODE 2947*	2,3,4,5,6	1.208,00

Armrests options

	CODE	DESCRIPTION	EUR €
❑	101	Clothing guard for armrest *Only i.c. with standard armrest CODE 106*	59,00
❑	4422	Armrest fixation in initial position *Only in combination with biomechanical armrest CODE 324 or swivelling armrest CODE 24*	63,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Armrest cushion LEFT

	CODE	DESCRIPTION	POSITION	EUR €
❑	975	Armrest padded (standard)	left	0,00
❑	4990	Armrest padded (standard) with swivel joint	left	392,00
❑	4982	Armrest hemiplegia lower arm support	left	440,00
❑	4993	Armrest hemiplegia lower arm support with swivel joint	left	468,00

Right armrest cushion

	CODE	DESCRIPTION	POSITION	EUR €
❑	974	Armrest padded (standard)	right	0,00
❑	4984	Armrest padded (standard) with swivel joint	right	392,00
❑	4981	Armrest hemiplegia lower arm support	right	440,00
❑	4991	Armrest hemiplegia lower arm support with swivel joint	right	468,00

Hand rests for hemiplegia forearm shell LEFT

	CODE	DESCRIPTION	POSITION	EUR €
❑	4988	Handrest flat for hemiplegia lower arm support	left	67,00
❑	4987	Handrest with ergonomic grip for hemiplegia lower arm support	left	67,00

⚠ *Not compatible with armrest padded on the left CODE 975*
Not compatible with armrest padded with swivel joint on the left CODE 4990

Hand rests for hemiplegia forearm shell RIGHT

	CODE	DESCRIPTION	POSITION	EUR €
❑	4985	Handrest flat for hemiplegia lower arm support	right	67,00
❑	4986	Handrest with ergonomic grip for hemiplegia lower arm support	right	67,00

⚠ *Not compatible with armrest padded on the right CODE 974*
Not compatible with armrest padded with swivel joint on the right CODE 4984

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Extension for hemiplegia lower arm support

	CODE	DESCRIPTION	POSITION	EUR €
❑	4989	Extension for hemiplegia forearm shell LEFT ⚠ *Not compatible with armrest padded on the left CODE 975* *Not compatible with armrest padded with swivel joint on the left CODE 4990*	left	89,00
❑	4983	Extension for hemiplegia forearm shell RIGHT ⚠ *Not compatible with armrest padded on the right CODE 974* *Not compatible with armrest padded with swivel joint on the right CODE 4984*	right	89,00

06 Leg supports

Side leg rests

	CODE	DESCRIPTION	EUR €
❑	293	80° Legrests, detachable and swing-away	0,00
❑	4092	Elevating leg rests (Manual), removable, Tool height-adjustable	690,00
❑	2086	Powered elevating, Legrests, removable,	2.390,00

Options for side leg supports

	CODE	DESCRIPTION	depth	LLL	EUR €
❑	808	Flip up footrests, angle adjustable			0,00
❑	57	One-piece footrest			452,00

Central legrest

	CODE	DESCRIPTION	LS-ANGLE	LLL	EUR €
❑	5092/5100	Central legrest, mechanically angle adjustable ⚠ *Only in conjunction with the footplate option for the central leg support, CODE 5109*	85°-95°	280-400 mm	640,00
❑	4935/5086	Electric central leg support (including calf pads, CODE 2945), electric length adjustment ⚠ *Only in conjunction with the footplate option for the central leg support, CODE 5109*	90°-170°	200-400mm	2.350,00
❑	4949/5086	Powered central legrest (incl. calf pad), electrically lowerable, length adjustment floor ⚠ *Only in conjunction with the footplate option for the central leg support, CODE 5109*	90°-170°	350-450mm	2.550,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Footplate options for central legrest

	CODE	DESCRIPTION	EUR €
❑	5109	Standard foot plate for central legrest	0,00

⚠ *Not compatible with lateral legrests*

Additional options Leg support

	CODE	DESCRIPTION	EUR €
❑	822	Toe strap (pair)	69,00
❑	823	Heel Loops (pair)	120,00
❑	2945	Calf pads for central leg support	0,00

07 Batteries/Chargers/Lighting

Batteries

	CODE	DESCRIPTION	EUR €
❑	5158	Batteries: 66 Ah (5H), 80 Ah (20H), maintenance-free	0,00
❑	700	Without batteries	0,00

Battery chargers

	CODE		EUR €
❑	732	Charger 8 A for 73 Ah (5H)	0,00
❑	733	Charger 12 A for 73 Ah(20h)	348,00
❑	701	Without charger	0,00

Lighting

	CODE	DESCRIPTION	EUR €
❑	932	Active LED-lighting front and rear acc. to StVZO	0,00
❑	703	Without active LED lighting	0,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

08 Control modules/mounting/brackets/adjustment modules

R-NET control system

Code	Control modules	120A power module, ISM8 adjustment module
702/415	R-NET LED	❑
702/416	R-NET LCD	❑ 990,00
702/398	R-NET LCD IR/BT	❑ 2.690,00
702/839/2859	R-NET OMNI II	❑ 3.017,00

⚠ *R-NET OMNI II: The selected control module is included*

Vehicle stability control module

	CODE	DESCRIPTION	EUR €
❑	918	Gyro module, vehicle stabilisation module ⚠ *Not compatible with the LED Code 415 control module*	1.490,00

Joystick Location

	CODE	DESCRIPTION	EUR €
❑	74	Assembly right	0,00
❑	60	Assembly left	0,00

Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Joystick /Control module holder

	CODE	DESCRIPTION	EUR €
❑	4859	Joystick bracket, length adjustable	0,00
❑	4860	Operating module holder, length and height adjustable, swingaway	210,00
❑	2074/4789	Attendant control mounted on the right, including 4-way bus distributor CODE 2859	793,00
❑	2060/4789	Attendant control mounted on the left, including 4-way bus distributor CODE 2859	793,00

External adjustment module

	CODE	DESCRIPTION	ASSEMBLY	EUR €
❑	4446/4597	External control panel TA, for 5 electrical functions ⚠ *Only in conjunction with control module bracket CODE 4860* *Not for use in conjunction with control module CODE 415*	Control module side	796,00
❑	4446/4599	External control panel TA, for 5 electrical functions ⚠ *Not compatible with control module CODE 415*	opposite the control module side	630,00

R-NET special control system

	CODE	DESCRIPTION	EUR €
❑	4660	Table top control can be swivelled to the right ⚠ *Not compatible with armrest CODE 24* *Not compatible with R-Net OMNI II CODE 839*	810,00
❑	4661	Table top control can be swivelled to the left ⚠ *Not compatible with armrest CODE 24* *Not compatible with R-Net OMNI II CODE 839*	810,00
❑	2238	Narrow rear headrest control, R-NET ⚠ *Only in conjunction with R-NET OMNI II CODE 839*	5.954,00
❑	2239	Wide rear head control, R-NET ⚠ *Only in conjunction with R-NET OMNI II CODE 839*	5.954,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Eye-Controlled driving system

	Item	DESCRIPTION	EUR €
❑	1096412	GazeDriver Set R-Net	34.998,00
❑	1096413	Therapy table - installation system	2.098,00
❑	1096410	Seat-installation system *Mandatory selection if GazeDriver is selected*	0,00

Operating module software

	CODE	DESCRIPTION	EUR €
❑	985	Series programming (standard)	0,00
❑	998	Special programming (deviation from standard)	0,00

External button

for users who can hardly operate the keys of the operating module or as an emergency/off function for accompanying person

	CODE	DESCRIPTION	EUR €
❑	596	Push-button, diameter 63 mm, yellow	259,00
❑	778	Push-button, diameter 30 mm, green	259,00
❑	779	Push-button, diameter 30 mm, yellow	259,00

Joystick attachments

	CODE	DESCRIPTION	CONTROL	EUR €
❑	553	TETRA fork vertical	R-NET	228,00
❑	4590	Joystick button ball handle 38 mm	R-NET	70,00
❑	4591	Joystick button T-handle	R-NET	146,00
❑	4593	Joystick button handle	R-NET	82,00
❑	4594	Joystick button handle	R-NET	74,00
❑	4595	Joystick knob softball 70 mm	R-NET	188,00
❑	4596	Golfball, small, 43 mm	R-NET	76,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

09 Accent colours

	CODE	FRAME	AKZENTE	EUR €
❑	186	Matt black	Magic Blue	0,00
❑	212	Matt black	Polar white	0,00
❑	227	Matt black	Red metallic	0,00
❑	180	Matt black	Lemon reflex	0,00
❑	214	Matt black	Sun yellow	0,00

10 Options around the wheelchair

	CODE	DESCRIPTION	EUR €
❑	4646	Therapy table swivelling sideways	696,00
❑	589	Shoulder Harness	378,00
❑	833	Safety belt (metal buckle)	59,00
❑	2859	Bus distributor 4-fold	83,00

11 Service

	Item	DESCRIPTION	EUR €
❑	55704	Flat fee for customised parameterisation (not eligible for discounts) ⚠ *The on-the-spot service covers only simple adjustment of parameter sets, such as acceleration, speed etc.*	556,00
❑	55709	Medically necessary parameterization to client-specific setting values ⚠ *The on-the-spot service covers simple program adjustments as well as the integration of special controls.*	1.124,00